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
## AGSS Forano

Dear customers,

Thank you for purchasing this Greggersen product.

If you have any questions or require information, please contact our sales and support team.

Tel.: +49-(0)40 739357-0  
E-mail: [sales@greggersen.de](mailto:sales@greggersen.de)

Please read the instructions for use thoroughly before using the medical device. Only use it once you have completely understood the application and function!

Medical device in  
accordance with:

DIN EN ISO 7396-2  
DIN EN ISO 9170-2  
MDR 2017-745

**Preface**

The AGSS Forano terminal unit from Greggerson Gasetechnik GmbH complies with the relevant provisions of Regulation (EU) 2017/745 of the European Parliament and Council on medical devices (MDR 2017/745) and is classified as a Class I active, non-invasive device.

**Intended use**

To transfer residual gases during inhalation anaesthesia at a suitable disposal point.

**Functionality**

Terminal units in anaesthetic gas scavenging systems (AGSS) are connected directly to a central compressed gas supply with a supply pressure of approx. 500 kPa compressed air or nitrogen. The ejector is connected to the gas supply and generates a negative pressure in the operational state. Excess anaesthetic gases are suctioned off by this negative pressure and together with the compressed gas are passed via the suction line into the outside air or into the exhaust air shaft of an air conditioning system.


The terminal unit is set to the operating state by connecting an AGSS connector or a hose with the appropriate connection. If the connection is disconnected again, the suction stops. An indicator is installed for operational control.

**Warnings**

- Prior to the use of the device, the user must be familiar with the product and the instructions for use and has to have understood them.
- The terminal unit must be checked for external damage or contamination before each use. It must be ensured that the terminal unit is correctly aligned and that the labelling on the terminal unit is horizontal.
- When cleaning, make sure that no liquids enter the terminal unit. After cleaning, the terminal unit must be dried completely.
- After installation of the terminal unit, an acceptance test according to DIN EN ISO 7396-2 must be performed.
- Only use the terminal unit with type 1L AGSS connections.
- All serious incidents must be reported to Greggerson and to the health authority responsible for the place of residence.

## Device description

Pos.	Description
1	O-ring 2,5x1,3
2	Closing cone AGSS
3	O-ring 9,25x1,78
4	O-ring 14x2
5	Washer 25x18x1
6	O-ring 16,5x1
7	Forano actuator AGSS
8	Connection carrier air
9	Connection indicator
10	Connection carrier air

### Accessories

- 902071 plug AGSS
- 902095 AGSS-plug, angled 45°

**Patient target group**

The patients themselves are not in direct contact with the product. The patient target group is therefore defined by the device connected.

When using the AGSS (anaesthetic gas scavenging system) Forano as intended, the focus is not on protecting the patient but rather on protecting persons in the immediate vicinity. The medical personnel in particular are affected.

**Intended users**

The medical device is intended for clinical use by medically qualified personnel (therapists, care personnel, doctors).

**Operation**

Before starting operation each time, an external visual inspection and a functional test have to be carried out. After visual inspection for obvious damage or contamination, the actuator must be activated and smooth movement verified.

To establish a connection, an AGSS plug or a hose with an AGSS type 1L connection is inserted into the terminal unit.

By pressing the actuator on the terminal unit, the connector is released from the terminal unit. The plug should be held by hand! The suction is interrupted.

**Cleaning**

The terminal unit can be wet cleaned with commercially available disinfectants and cleaning agents. When cleaning, make sure that no liquid enters the terminal unit. The terminal unit must be completely dry prior to subsequent use. The hygiene plan of the healthcare facility must be observed.

**Installation**

- The installation housing in the hollow-wall / under plaster must first be secured in the wall!
- Hollow-wall installation: Fasten the housing with six screws left and right, solder the pipes to the piping system
- Under plaster: Plaster the housing in the wall, solder the pipes to the piping system
- Connect the terminal unit to the pipe or ejector connections using the connecting lines (ejector can be plugged and fixed with a screw!)
- Push the terminal unit with the front panel into the housing and fasten it with two screws
- Supply systems: Fix the terminal unit with two screws, insert the ejector into the exhaust air line and secure it with a nut
- Connect the terminal unit, ejector and pipe connections to the connection lines
- Adjust suction capacity with valve spindle (11 o'clock position) with hexagon fittings (2,5 mm)
- Press the actuator on the terminal unit
- Check anesthesia gas exhaust system according to DIN EN ISO 7396-2

### Maintenance and inspection

Always perform a visual inspection for external damage before using the terminal units (external damage). A function and tightness check should be performed at least once a year (if necessary: replacement of wearing parts), maintenance should be carried out every six years (measurement and safety check) in particularly high-use areas more frequently!

Faults and damage must only be dealt with by authorised expert personnel. Only use genuine GREGGERSEN replacement parts when performing repairs.

### Disposal

The terminal unit is to be disposed of in accordance with the relevant national regulations.


### Technical Data

Design:	All-metal design, pluggable ejector for vacuum generation, single-handed operation for connecting/disconnecting, indicator for operational control, integrated valve for throughput adjustment.
Connection:	<u>Supply lines underplaster/ hollow wall/ onplaster</u> Pressure gas line: Cu pipe 8x1 Exhaust line: Cu pipe 15x1 <u>Terminal unit / ejector</u> Compressed gas supply: Connection for 4 mm PA hose or 7 mm barbed hose fitting Ejector: Connection for 4 mm PA hose or 7 mm barbed hose fitting Anaesthetic gas dissipation: 15 mm spiral hose Indicator: Connection for 4 mm PA hose
Throughput:	min. 50 l/min at 500 kPa compressed air
Gas consumption:	approx. 12 l/min at suction throughput of 50 l/min
Installation:	Concealed / hollow-wall installation Ceiling supply systems Wall-mounted supply systems


Our products are partly made of brass with a lead content of more than 0.1%. As an alloy component, lead is firmly bound in the material so that no exposure is to be expected.


### Installation dimensions hollow wall/ under plaster


### Installation dimensions on-plaster


### Installation dimensions ceiling support bed head units

### Installation dimensions ceiling pendants

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-29 / 21031 Hamburg / Germany