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

Operating alarm unit

Dear customers,
thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggersen.de

Usage

Monitoring and display of the pressure status of the central gas supply systems. With alarm message in the event of impermissible pressures.

Structure

The front is coated white (RAL 9016). The control panel is gray (RAL 7040). In the housing of the operating signal there are up to 6 monitoring channels. Each channel of the M2 module is equipped with 3 LEDs (red, green and red), each channel of the M12 module with 2 LEDs (yellow and green). The operating signal is connected to the peripheral device via their sensor connection (2-wire cable).

Operation

Acknowledge alarms The active alarm can be silenced by pressing the „Reset“ key. This means that the acoustic alarm is switched off for a period of 12 minutes. The displays / relays remain in the alarm state. If the cause of the alarm is still present after 12 minutes have elapsed, the audible alarm is reactivated. If multiple alarms are active at the same time in a multiple alarm situation, e.g. Gas1 too high pressure, Gas2 too low pressure, then pressing „Reset“ will silence the audible alarm for all active alarms. A newly added alarm immediately activates the acoustic alarm again.

Technical data

Type:	Display of up to six alarm messages with three LEDs each (red, green, red), or 2 LEDs each (yellow, green), potential-free contact (changeover contact) per message. Surface-mounted housing made of sheet steel painted white (RAL 9016), flush-mounted housing frame made of sheet steel painted white (RAL 9016), built-in housing made of stainless steel. Control panel gray (RAL 7040). Electrical connection top or bottom, integrated power supply.
Dimensions:	Underplaster / Hollowall Housing (HxWxT) ca.: 123 x 238 x 65 [mm] Cover (HxWxT) ca.: 131 x 303 x 16 [mm] Surface-mounted housing (HxWxT) ca.: 131 x 303 x 65 [mm]
Power supply unit: (integrated)	Input: 100...240 V AC, 50-60 Hz Output: 24V DC
Elektronics:	Input voltage: 12...24V AC / DC
Relay contact:	Changer, 24V DC 1A
Guideline:	DIN EN ISO 7396-1 DIN EN ISO 60601-1-8

Important notes

- Connect the ground connections properly.
- The power supply unit may only be connected by a qualified electrician.
- Observe the enclosed documentation of the power supply unit.
- There are no warranty claims for damage or malfunctions caused by improper transport, improper storage or improper use.

Recommended installation sequence

Mounting

- Mounting of the rear part in/on the wall.
- Connect the sensors to the electronics
- Connect the shut-off box to the supply mains.- Connection of the power supply unit to the mains supply must be carried out by a qualified electrician.
- Connect the cable of the power supply unit to the warning.
- For flush mounting/cavity wall mounting: Screw the cover frame back onto the rear part using the four M4 screws, aligning horizontally.
- Insert the carrier plate with electronics into the frame (first at the top then push it downwards, see pictures on page 10).- Pay attention to the cables (do not forget the grounding cable).
- Staple the magnetic foil to the carrier plate.
- Switch on the power supply for the power pack.

Testing

- Check the labeling of the signals
- Complete function check

Maintenance / Service

Every Operating alarm unit must be subjected to a regular visual inspection (external damage). A function test must be carried out at least once a year. To do this, press the test button (for the procedure, see page 9). Malfunctions and damage may only be remedied by authorized specialist personnel. Only use original GREGGERSEN spare parts for repairs!

Cleaning


Rub the front with a slightly damp, soft cloth. Make sure that no cleaning fluid runs into the housing.

Attention: Under no circumstances use cleaning agents with abrasive additives or contaminated cleaning water with granular components !!!

Terminal assignment auxiliary warning

Pos.	Description
1	Relay contact
2	Sensor connection: 1→24V DC (white); 2→Signal (brown)
3	Supply voltage
4	Free
5	Free
6	Interface programming RS485
7	Relays contact main board fault and door alarm
8	Free
9	Connection "M8" board (Option)

current condition emergency alarm unit	relay contact (NI - NC) (dry contact)
currentless	open
normal operating (green)	close
alarm (red)	open

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Error messages

No.	Description	Prio	LED red	LED green	LED red	Signal
1	normal operating	0	off	on	off	off
2	test alarm	1	flashing	off	flashing	on
3	pressure to high	10	off	off	flashing	on
4	pressure to low	10	flashing	off	off	on
5	door opened	2	off	flashing	off	on
6	pressure switch alarm	10	flashing	off	flashing	on
7	error in module M1	9	flashing	flashing	flashing	off
8	error in module M2/M3	9	off	off	off	on
9	CRC error in EEPROM M1	7	off	flashing	on	on
10	CRC Error in EEPROM M2/M3	7	on	flashing	off	on
11	communication error M1->M2/M3	8	off	off	off	on
12	communication error M2/M3->M1	8	flashing	flashing	flashing	off
13	communication error separate alarm unit	3	on	flashing	on	off
14	short circuit at sensor M2/M3	9	on	flashing	flashing	on
15	Sensor cable break on M2/M3	9	flashing	flashing	an	on
16	loss communication main alarm unit -> separate alarm unit	5	on	flashing	on	on
17	Aloss communication with M8	4	blinken	on	on	on
18	error M8	8	on	on	on	off

Test alarm sequence

By pressing the „Test“ button, the following sequence self-test of the EEPROM and program (FLASH) memory of the microcontroller of the operating signal begins

1. All LEDs on for two seconds;
2. All LEDs off for one second;
3. Switching on the acoustic signal generator; With M2 module: All red LEDs start to flash; With M12 module: All yellow and green LEDs start to flash; The self-test switches off after 10 seconds.

Any real alarm immediately interrupts the test sequence.

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