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
Instructions for use

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Automatic switch cabinet

HU10

Deutsch.....2
English.....7

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Dear customers,

Thank you for purchasing this Greggersen product.

If you have any questions or require information, please contact our sales and support team.


Tel.: +49-(0)40 739357-0
e-mail: sales@greggersen.de

1. Warnings

- Prior to the use of the device, the user must be familiar with the product and the instructions for use and has to have understood them.
- Only trained staff may work on the system.
- Never enter the cylinder bank room with a naked light or flame.
- It is essential that oil and grease are avoided in the control centre! **EXPLOSION HAZARD!**
- Secure all cylinders with a chain to protect against their falling over.
- Carry out a once-daily inspection of cylinder banks in continuous operation.
- Use only original spare parts from the manufacturer.
- Open the cylinder valves slowly, so as to avoid pressure surges.
- UVV (Accident Prevention Regulations) oxygen must be observed when using the system.
- The system meets DIN 13260

2. Structure

DM1	HP-pressure reducer, right side
DM2	HP-pressure reducer, left side
DM3	Pressure reducer, working reducer
M1	Pressure gauge, right cylinder bank
M2	Pressure gauge, left cylinder bank
M3	Pressure gauge, mean pressure
M4	Pressure gauge, working pressure
SV	Safety valve, 7 bar
KN	Contact switch low pressure
V1	Main shut-off valve
H1	Lever
EL	Input left
ER	Input right

3. Description

The Automatic switching system is a component of the central gas supply system for technical gases. With two-way cylinder bank supply, they are designed for capacities of up to 10 m³/h.

In principle, the switching system has the following tasks:

- To reduce the pressure in the cylinder banks to a constant working pressure.
- To switch to the other cylinder bank if the pressure drops below a specific level.
- To ensure the replacement of an empty cylinder bank without the gas supply being interrupted.
- To record the correct pressure ratios in the system.

All the equipment for pressure reduction, switching and monitoring are accommodated in a splash-proof housing. The pressure gauges in the central section of the system display the pressure in the right (M₁) and left (M₂) cylinder banks and the working pressure (M₄).

The gas is passed to the switching system through manifolds from the cylinder banks. The cylinder pressure is then reduced in two stages to a working pressure which is normally 5 bar. The safety valve SV opens if the working pressure exceeds 7 bar.

4. Function

The pressure reducers M₁ and DM₂ are supplied with gas from the relevant cylinder bank and reduce pressure to a mean of about 6 bar. The eccentric control lever is used to select that side which is to start up (the lever must point in the direction of the pressure reducer to be triggered). This causes the outlet pressure of the corresponding pressure reducer to increase by 2 bar, the pressure difference closes the backflow check valve on the other side.

The lever at DM₃ can be used to set the desired working pressure (normally 5 bar) and this pressure can be read off at M₄. A safety valve opens at 7 bar. When the stop valve at the outlet is open, the piping network is supplied with gas.

5. Operation

In normal operation, the switching system is supplied by one cylinder bank while the other side is available as a reserve. The lever on the eccentric normally points in the direction of that cylinder bank which is operating. M₁ and M₂ show the cylinder pressure, M₄ the working pressure.

If the cylinder pressure drops to 6 bar, the system automatically switches to the other side. Before the empty side can be replaced, the eccentric control lever must be turned 180° (pointing in the direction of the other side providing the gas supply) since otherwise the new cylinders will immediately take over the supply of gas again.

After the connection between the empty bank and the switching system has been broken with the relevant stop valve, the cylinders can be replaced and the stop valve reopened.

6. Recommended procedure for replacing cylinders

- Close cylinder valve and main stop valve.
- Vent the line with the vent valve where manual connection elbows are used.
- Undo connection nut on cylinder valve.
- Fit protective cover on cylinder and replace with full cylinder.
- Remove protective cap from cylinder and screw connection nut on to cylinder valve (screw on without tools for manual connection nuts).
- Turn the eccentric control lever in the switching system.
- Open the cylinder valve slowly and then open the main stop valve slowly.


7. Maintenance

The entire switching system should be visual inspected by trained specialist staff at regular intervals. This inspection includes, among other things, a check of the correct pressure ratios (M₁ to M₄) and the tightness of the pipe connections.

Maintenance of the system with replacement of all wear parts should be carried out every 5 years by authorized specialist personnel. Only use original **GREGGERSEN** spare parts for repairs and maintenance.

7. Spare parts

325265	Repair set for switch cabinet (HU), 1st step
325266	Repair set for switch cabinet (HU), 2nd Step
900466	Repair set safety valve 7 bar
900467	Repair set safety valve 11 bar

www.greggersen.com

GREGGERSEN Gasetechnik GmbH / Bodestr. 27-29 / 21031 Hamburg / Germany

23.06.2021/M.Oels/ 104274-01 DOKU - HU10 de, en