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INSTRUCTIONS FOR USE

Laborkleinanlage Small Lab Unit

Deutsch	2
English	4

Laborkleinanlage

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Small Lab Unit

Dear customers,

Thank you for purchasing our Greggersen product.
If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739357-0, sales@greggersen.de

Please read the instructions for use thoroughly before using the medical device. Only use it once you have completely understood the application and function!
All serious incidents must be reported to Greggersen Gasetechnik GmbH.

Symbols:

LOT Date of manufacture
calendar week year

Manufacturer

Keep free from
oil and grease

Usage

Reduces cylinder pressure to a constant working pressure.

Components

The set includes all necessary components for easy and safe installation:

- Mounting bracket for stable fixation
- Connection elbow for simple connection
- Shut-off valve for safe control of gas flow
- Pressure regulator with content and working pressure gauges for precise pressure monitoring
- Fittings for tight and reliable connections

Installation and Commissioning

1. **Mount the high-pressure valve**
Secure the high-pressure valve to the wall using the mounting bracket. The recommended mounting height is 1.70 m, but it can be adjusted based on the size of the gas cylinder.
2. **Connect the high-pressure elbow**
Firmly screw the high-pressure elbow to the high-pressure corner valve. Securely connect the other side to the gas cylinder. Ensure cleanliness and proper sealing of all connections. For oxygen applications, always work oil-free!
3. **Install the pressure regulator**
Attach the pressure regulator to the high-pressure corner valve. Ensure the connection is tight and secure.
4. **Prepare the pressure regulator**
Before commissioning, fully release the adjustment screw on the pressure regulator until no spring pressure is felt.
5. **Open valves and check pressure**
Slowly open the cylinder valve and the high-pressure corner valve. The content pressure gauge on the pressure regulator will display the cylinder pressure. At this point, gas does not yet flow. Check the entire system for leaks, e.g., by brushing connections with soapy water or observing for pressure drops after closing the shut-off valve.
6. **Set the working pressure**
Once no leaks are detected, adjust the desired working pressure by turning the adjustment screw on the pressure regulator. Then, open the shut-off valve downstream of the pressure regulator.

Decommissioning

Decommission the system in reverse order:

1. Close the shut-off valve downstream of the pressure regulator.
2. Close the high-pressure corner valve.
3. Close the cylinder valve.

Cylinder Replacement

1. Close the cylinder valve and the high-pressure corner valve.
2. Loosen the nut of the connection elbow at the cylinder valve.
3. Replace the cylinder and reconnect as described above.

Always follow applicable safety regulations and work carefully, especially when handling oxygen.

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany