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Switch cabinet MediControl MC2025P

Dear customers,

thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggersen.de

Definitions

Main supply / reserve supply

- In a distribution system with cylinder batteries, the two cylinder batteries form the main supply source and another cylinder battery acts as reserve supply.

Use

The central gas distribution system **MediControl** guarantees a continuous supply with medical or technical gases (oxygen, nitrous oxide, carbon dioxide and other gases). In line with DIN EN ISO 7396-1, the system is rated for 3 supply sources. It consists basically of 3 components: the pneumatic switch-over system with integrated contact-switch component, a reducer panel and the manifolds.

The pneumatic switch-over system controls and monitors the supply of the gas distribution system. At the same time it reduces the pressure of the cylinder batteries.

The reducer panel reduces the pressure of the reserve supply (for 2 cylinder batteries + 1 reserve supply) on the network pressure.

The system can be configured in two different versions:

- 2 cylinder batteries + 1 reserve supply
- 2 cylinder batteries as reserve supply

Important notes

- The person using/operating the unit has to ensure that any operators are familiar with the unit's operating instructions and function.
- There is no entitlement to any warranty claims in the event of damage or malfunction arising from improper transport, improper storage or other than the intended use.
- This unit must not be modified without the manufacturer's permission!
- Only use oxygen applicable lubricants (fire/explosion hazard)!
- Any changes to the setting of the safety valve will change the operating pressure and the accessories can be damaged.
- Always open valves slowly! Any sudden opening may cause a fire and explosion hazard; surges in pressure can damage the unit or its accessories.


- At the end of its service life: arrange for proper disposal of the unit following consultation with the disposal company responsible. Observe the applicable legal regulations.
- Also include the attachment of the copper piping's earthing connections (clamps) into the earthing concept.

Description of functions

The central gas distribution system MediControl is intended to safeguard the continuous supply with medical or technical gases.

The gas is brought from the cylinder batteries through the manifolds to the switch-over system. The high pressure on the cylinder batteries is reduced in two stages in the switch-over system. The two high-pressure reducers (DM1, DM2) reduce the high pressure from the cylinder batteries to a medium pressure. The left-hand high-pressure reducer (DM1) is adjusted to a pressure of approx. 850 kPa, the right-hand one to approx. 750 kPa. This guarantees priority for the reducer on the left side. Each high-pressure reducer is equipped with a safety valve (SV1, SV2) which opens at a pressure of approx. 1100 kPa.

left side priority

right side priority

The medium pressure is reduced by two low-pressure reducers (DM3, DM4) to the network pressure. Each reducer has an output which corresponds to the rated output, i.e. even when one of the reducers is removed or being serviced, the switch-over system still performs to the rated output level. When the network pressure exceeds 700 kPa, a safety valve opens (SV3).

An emergency supply point (NE1) allows for the supply network to be supplied via a NIST-connector coupling in an emergency. Gas supply can then be provided for example by means of a gas cylinder. Please note that the corresponding shut-off valves of the system are closed during an emergency supply (Control lever in the direction of the flow: valve open), and that after an emergency supply, the system has to be started up again according to the chapter „initial commissioning“.

The lever (which is mount on one high-pressure-reducer) is used to select that side which is to start up (the lever must point in the direction of the pressure reducer to be triggered!). This causes the outlet pressure of the corresponding pressure reducer to increase by 2 bar, the pressure difference closes the one-way valve on the other side.

Technical data

pneumatic switch-over system

Design:	Twin-designed, two-stage network pressure reduction with pneumatic priority switching, removable painted steel plate housing with easily opened hood, safety valves for medium pressure and network pressure, emergency supply point (NIST)
Operating pressure:	Input max. 20.000 kPa (200 bar) Output 450...550 kPa (4,5...5,5 bar), max. 800 kPa (max. 8 bar)
Dimensions:	Measurements (WxHxD) approx. 380 x 840 x 300 [mm]
Weight:	approx. 35 kg
Capacity:	25 Nm ³ /h
Incoming:	G3/4" a flat (cylinder batteries)
Outgoing:	to the distribution panel: copper pipe Ø 22 mm Safety valves: copper pipe Ø 12 mm
Guidelines:	DIN EN ISO 7396-1 DIN EN ISO 10524-2 DIN VDE 0107 Medical Device Directive 93/42/EEC


contact-switch

Design:	oil and grease-free
Incoming pressure:	DS1, DS2: 20.000 kPa (200 bar) DS3: 1.000 kPa (10 bar)
Voltage max:	DS1, DS2: 250 volt DS3: 250 volt
Power max:	DS1, DS2: 5 amp DS3: 5 amp

Schematic structure of the switch-over system

grid

Key for the schematic structure

Abb.	Description
DM ₁ /DM ₂	High-pressure reducer
DM ₃ /DM ₄	Low-pressure reducer
DMBGR ₁ /DMBGR ₂	Reducer assembly
DS ₁ /DS ₂	High pressure contact switch cylinder manifold
DS ₃	2 pressure switch network-pressure (1 normally open switch, 1 normally closed switch)
F ₁ /F ₂ /F ₃ /F ₄	Filter
M ₁ /M ₂	Pressure gauge 0-31.500 kPa
M ₃ /M ₄	Pressure gauge 0-1.600 kPa
NE ₁	Emergency supply point (NIST)
RV ₁ /RV ₂	One-way valve medium-pressure
RV ₃ /RV ₄	One-way valve low-pressure
V ₁ -SV ₂	Safety valve medium-pressure, 1.100 kPa
SV ₃	Safety valve network-pressure, 700 kPa
V ₁ -V ₂	Shut-off valves, low-pressure
V ₃	Main shut-off valve cylinder supply

Installation instructions

Space required

A distribution system with gas cylinders should be installed in a specially designed room, or a room which has been suitable modified, which is well ventilated and fire-proof.

Recommended installation procedure


- Attachment of the switch-over system to the wall.
- Attachment of the manifolds to the wall.
- Connection of the manifolds to the switch-over system.
- Route the vent pipes of the safety valves and the manifolds outside:
->The vent pipes for the medium pressure stage and those for the network pressure **must** be routed separately.
- After installation, rinse the whole system.
- Route and connect up electrical connections:
-> Connect the clamps „Common“ and „Norm Closed“ of the contact switch for the left manifold (DS₁) with the clamps 1 and 8 (plug 6) of the Greggersen alarm unit. See also documentation alarm units.
-> Connect the clamps „Common“ and „Norm Closed“ of the contact switch for the right manifold (DS₂) with the clamps 1 and 7 (plug 6) of the Greggersen alarm unit. See also documentation alarm units.
-> The both contact switches for network-pressure(DS₂) must be connect in series (connect one clamp of the normally open switch with one clamp of the normally closed switch). The other both clamps of the contact switches connect at the clamps 1 and 9 (plug 6) of the Greggersen alarm unit. See also documentation alarm units.
- **!!! Caution !!!** During soldering work, the pipes **must** be rinsed with inert gas!

Installation instructions

- Installation may only be carried out by authorised qualified staff.
- Only copper pipes as per DIN EN 13348 may be used for installation.
- Reports must be kept of the acceptance and handing over procedures for the gas distribution system.
- All pipes and connections must always be kept free of oil and grease! **!!! Risk of explosion !!!**
- To connect the cylinders with the high pressure header use a tool to hold the high pressure check-valve! New check-valves shall be installed with a torsion force of appr. 40 Nm.

Structure of the central gas distribution system

Initial commissioning

- Before initial commissioning, the system must be vented and **all** shut-off valves must be **closed**.
- Installation must be finished completely.

1. Slowly open the main shut-off valves of the manifolds

- Observe the pressure at pressure gauges M₃, it must remain constant and adjust to a pressure of approx. 850 kPa .
- The pressure on the left-hand side should be approx. 200 kPa higher than on the right-hand side. If the pressure increases slowly, the reducer must be checked.

2. Open shut-off valve V₁ to V₂

- Observe the pressure at pressure gauge M₄, it must adjust to the network pressure and remain constant.
- If the pressure increases slowly, the reducer must be checked.

3. Slowly open shut-off valve V₃

- Observe pressure at pressure gauge M₄, it must adjust to the network pressure.

Maintenance / Inspection Inspection

The whole gas distribution system should be subject to visual inspection at regular intervals by an authorised qualified technician, checking the functions and pressure conditions in the system.

Maintenance

Maintenance of the system is recommended once every twelve months and may only be carried out by an authorised qualified technician. This includes among others, checking the correct pressure conditions and tightness of the pipe connections together with the valves. All service parts should be replaced every 5 years.

If individual components have to be removed for maintenance or repair, the corresponding shut-off valves before and after the component have to be closed. The redundant design with double components in the switch-over system means that the network continues to be supplied with gas, so that it is not necessary to adjust the reducers.


If the system has to be decommissioned for a brief time (major repairs, replacement of whole components), it is possible for gas to be supplied by the emergency supply. Please make sure that the main shut-off valve of the system is closed during emergency gas supply. To start up the system again for normal operation, please proceed according to the chapter "Initial commissioning". Malfunctions and damage may only be rectified by authorised trained personnel. Use only original GREGGERSEN service-parts for repairs!

Service-Parts

No.	Description
325313	Maintenance kit MC2025P „all in one“
325284	Maintenance kit for low-pressure reducer MC2025P
325285	Maintenance kit for high-pressure reducer MC2025P
325277	Maintenance kit for connections MC2025P
800660	CU seal high-pressure header G3/4"

Explosion drawing

Legend explosion drawing

Abb.	Assembly
AN	Outgoing to the network AVL Shut-off valve left
AVR	Shut-off valve right
DM1	Left high-pressure reducer with one-way valve
DM2	Right high-pressure reducer with one-way valve
DM3	Left low-pressure reducer
DM4	Right low-pressure reducer
DS1	Pressure contact switch left cylinder battery
DS2	Pressure contact switch right cylinder battery
DS3	Pressure contact switch network pressure
ML	Pressure gauge row
EL	Incoming, left cylinder manifold
ER	Incoming, right cylinder manifold
HV	Main shut-off valve
M1	Left high pressure gauge
M2	Right high pressure gauge
M3	Medium pressure gauge
M4	Network pressure gauge
NE1	Emergency supply point (NIST)
SVMD	Safety valve medium pressure, 1.100 kPa
SVND	Safety valve network pressure, 700 kPa
ZMD	Junction medium pressure

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