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Switch Cabinet MediControl MC2050 / MC2100

Dear customers,

thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggersen.de

Definitions

Main supply / reserve supply

- In a distribution system with tank, the tank is the main supply source and the cylinder batteries (bundles of cylinders) form the reserve supply (emergency supply).
- In a distribution system with cylinder batteries, the two cylinder batteries form the main supply source and another cylinder battery acts as reserve supply.

Reducer panel

- In a distribution system with tank, the reducer panel reduces the pressure of the tank on the network pressure.
- In a distribution system with cylinder batteries, the reducer panel reduces the pressure of the reserve supply on the network pressure.

Use

The central gas distribution system **MediControl** guarantees a continuous supply with medical or technical gases (oxygen, nitrous oxide, carbon dioxide and other gases). In line with DIN EN ISO 7396-1, the system is rated for 3 supply sources. It consists basically of 3 components: the electronic switch-over system with integrated electronic component, a reducer panel and the manifolds. The electronic switch-over system controls and monitors the supply of the gas distribution system. At the same time it reduces the pressure of the cylinder batteries. The status of the whole system can be checked at any time on a graphic display. The reducer panel reduces the pressure of the reserve supply (for 2 cylinder batteries + 1 reserve supply) on the net-work pressure.

The system can be configured in three different versions:

Tank + reserve supply (2 cylinder batteries)

- 2 cylinder batteries + 1 reserve supply
- 2 cylinder batteries as reserve supply

Important notes

- The person using/operating the unit has to ensure that any operators are familiar with the unit's operating instructions and function.
- There is no entitlement to any warranty claims in the event of damage or malfunction arising from improper transport, improper storage or other than the intended use.

- This unit must not be modified without the manufacturer's permission!
- Only use oxygen applicable lubricants (fire/explosion hazard)!
- Any changes to the setting of the safety valve will change the operating pressure and the accessories can be damaged.
- Always open valves slowly! Any sudden opening may cause a fire and explosion hazard; surges in pressure can damage the unit or its accessories.
- At the end of its service life: arrange for proper disposal of the unit following consultation with the disposal company responsible. Observe the applicable legal regulations.
- Also include the attachment of the copper piping's earthing connections (clamps) into the earthing concept.

Description of functions

The central gas distribution system **MediControl** is intended to safeguard the continuous supply with medical or technical gases.

The gas is brought from the cylinder batteries through the manifolds to the switch-over system. The high pressure on the cylinder batteries is reduced in two stages in the switch-over system. The two high reducers (DM1, DM2) reduce the high pressure from the cylinder batteries to a medium pressure. The left-hand high reducer (DM1) is adjusted to a pressure of approx. 850 kPa, the right-hand one to approx. 750 kPa. This guarantees priority for the reducer on the left side, even in the event of a power failure. Each high reducer is equipped with a safety valve (SV1, SV2) which opens at a pressure of approx. 1100 kPa.

The medium pressure is reduced by two low reducers (DM3, DM4) to the network pressure. Each reducer has an output which corresponds to the rated output, i.e. even when one of the reducers is removed or being serviced, the switch-over system still performs to the rated output level. When the network pressure exceeds 600 kPa, a safety valve opens (SV3).

An emergency supply point (NE1) allows for the supply network to be supplied via a NIST-connector coupling in an emergency. Gas supply can then be provided for example by means of a gas cylinder. Please note that the corresponding shut-off valves of the system are closed during an emergency supply (Control lever in the direction of the flow: valve open), and that after an emergency supply, the system has to be started up again according to the chapter „initial commissioning“.

The electronic switch-over system **Aeolus** controls and monitors the central gas distribution system. The current status of the system is shown on a display. Pressure sensors measure the pressure of the cylinder batteries (DS1, DS2) and the network pressure (DS3). Two pneumatic valves (PV1, PV2) are responsible for controlling the gas flow of the two cylinder batteries.


There are pressure gauges at all important control points so that the state of the system can still be assessed even in the event of a power failure.

Technical Data

Design:	Twin-designed, two-stage network pressure reduction with pneumatic priority switching, removable painted steel plate housing with easily opened hood, safety valves for medium pressure and network pressure, emergency supply point (NIST).
Operating pressure:	Input max. 20000 kPa (200 bar) output 450...550 kPa (4,5...5,5 bar), max. 800 kPa (max. 8 bar)
Dimensions:	(WxHxD) ca. 480 x 1100 x 330 [mm]
Weight:	ca. 45 kg
Capacity:	50 Nm ³ /h or 100 Nm ³ /h depending on the version
Incoming:	G3/4" a flat (cylinder manifolds)
Outgoing:	Distribution: copper pipe Ø 22 mm
Safety valves:	copper pipe Ø 12 mm
Guidelines:	DIN EN ISO 7396-1 DIN EN ISO 10524-2 DIN VDE 0107 Medical Device Directive 93/42/EWG
Elektronik:	Multiple line graphic status display, programmable, galvanized metal housing, terminal strip on electronic housing, sensor technique, Logbook, external power supply
Power supply:	Aeolus: 24 DC or 24 AC power adapter: input: 100 ... 240 V AC, 50 ... 60 Hz power adapter: output: 24 V DC, 1,75 A power consumption: < 30 W
Sensors:	grid / tank: 0...1600 kPa (0...16 bar)
High-Pressure:	0...25000 kPa (0...250 bar)
Supply voltage:	10...30 V DC
Signal:	4...20 mA
Outgoing:	9 potential-free contacts, max. 1 A / 60V
Magnet valves:	24 DC, 4W

Schematic structure of the switch-over system

Key for the schematic structure

Short cut	Description
DM ₁ /DM ₂	High-pressure reducer
DM ₃ /DM ₄	Low-pressure reducer
DMBGR ₁ /DMBGR ₂	Reducer assembly
DMT ₁	Reducer panel
DS ₁ /DS ₂	Pressure sensor cylinder battery, 0-25000 kPa
DS ₃	Pressure sensor network pressure, 0-1600 kPa
F ₁ /F ₂ /F ₃ /F ₄	Filter
M ₁ /M ₃	Pressure gauge 0-31500 kPa
M ₂ /M ₄ /M ₅	Pressure gauge 0-1600 kPa
NE ₁	Emergency supply point (NIST)
PV ₁ /PV ₂	Pneumatic valve
SV ₁ -SV ₂	Safety valve medium pressure, 1100 kPa
SV ₃	Safety valve network pressure, 600 kPa
V ₁ -V ₄	Shut-off valves, medium pressure
V ₅ -V ₆	Shut-off valves, low pressure
V ₇	Main shut-off valve cylinder supply
V ₈	Main shut-off valve cylinder supply

Installation instructions

Space required

A distribution system with gas cylinders should be installed in a specially designed room, or a room which has been suitable modified, which is well ventilated and fire-proof.

Recommended installation procedure

- Attachment of the switch-over system to the wall
- Attachment of the manifolds to the wall
- Attachment of the reducer panel
- Connection of the manifolds to the switch-over system
- Connection of the reducer panel to the switch-over system
- Route the vent pipes of the safety valves and the manifolds out-side
The vent pipes for the medium pressure stage and those for the network pressure **have to** be routed separately
- After installation, rinse the whole system
- Install the casing with the power supply unit on the wall next to the plant
- Route and connect up the electrical connections
Connect the pressure sensor of the reducer panel to the electronic component
Follow the documentation of the external power supply
!!! Connecting while the power is switched off !!!
- **Slowly** open the valves of the manifold and check the pressure
- Suspend the housing in the mounting plate

- Turn the power supply on

! Electrical safety test

The electronics may only be installed and commissioned by qualified technicians (electricians). Check the product according to DIN EN 62353 prior to initial commissioning and for every re-commissioning. Always follow state-specific regulations.

Check whether the PE conductor connection is properly and securely attached.

Measuring the PE conductor resistance:

1st measurement: electronics to ground.

2nd measurement: mounting plate to ground.


Do not exceed the threshold ≤ 0.3 ohms (according to DIN EN 62353, VDE 0751:2008).

- Function check for all components
- **!!! Caution !!!** During soldering work, the pipes must be rinsed with inert gas

Installation notes

- Installation may only be carried out by authorised qualified staff.
- Only copper pipes as per DIN EN 13348 may be used for installation.
- Reports must be kept of the acceptance and handing over procedures for the gas distribution system.
- All pipes and connections must always be kept free of oil and grease! **!!! Risk of explosion !!!**
- To connect the cylinders with the high pressure header use a tool to hold the high pressure check-valve! New check-valves shall be installed with a torsion force of approx. 40 Nm.

Structure of the gas supply system

Initial commissioning

- Before initial commissioning, the system must be vented and **all** shut-off valves must be **closed**
- Installation must be finished completely
- The electronic component must **not** be in operation

1. Slowly open the main shut-off valves of the manifolds

- Observe the pressure at pressure gauge M2 and M4, it has to remain constant and adjust to a pressure of M2 = approx. 850 kPa and M4 = approx. 750 kPa.
- The pressure on the left-hand side should be approx. 100 kPa higher than on the right-hand side. In the event of inadequate pressure difference, it is possible to feed both cylinder batteries into the supply network at the same time with the electronic component switched off.
- If the pressure increases slowly, the reducer must be checked.

2. Open shut-off valve V1 to V6

- Observe the pressure at pressure gauge M5, it must adjust to the network pressure and remain constant.
- If the pressure increases slowly, the reducer must be checked.

Operation with 3 supply sources

3. Slowly open shut-off valve V9

- Observe pressure at pressure gauge M6, it must adjust to the network pressure
- The pressure at pressure gauge M6 should be approx. 50 kPa higher than at M5
-

4. Slowly open the main shut-off valve V8

so that the pipe network is filled with gas

- Observe the pressure at pressure gauge M6 until the whole pipe network has filled and the pressure has adjusted to the network pressure
-

5. Open shut-off valve V7

- in distribution systems with a tank, it is possible that when the electronic component is switched off, the left hand side of the system also works too, i.e. the left-hand cylinder battery also supplies gas into the supply network.

Operation with 2 supply sources

3. Slowly open main shut-off valve V7

- Observe pressure at pressure gauge M5, until the pipe network has filled and the pressure has adjusted to the network pressure.

6. Switch on the power supply for the electronic Aeolus.

- No fault message may appear on the display if the electronic component.

Inspection

The whole gas distribution system should be subject to visual inspection at regular intervals by an authorised qualified technician, checking the functions and pressure conditions in the system.

Maintenance / Inspection

Maintenance of the system is recommended once every twelve months and may only be carried out by an authorised qualified technician. This includes among others, checking the correct pressure conditions and tightness of the pipe connections together with the valves. All service parts should be replaced every 5 years. If individual components have to be removed for maintenance or repair, the corresponding shut-off valves before and after the component have to be closed. The redundant design with double components in the switch-over system means that the network continues to be supplied with gas, so that it is not necessary to adjust the reducers. For maintenance of the reducer panels, first of all the shut-off valve of the tank (V9) has to be closed. This means that the left hand or right-hand cylinder battery then takes over the gas supply. It is then possible to close the main shut-off valve (V8). If the system has to be decommissioned for a brief time (major repairs, replacement of whole components), it is possible for gas to be supplied by the emergency supply. Please make sure that the main shut-off valve of the system is closed during emergency gas supply. To start up the system again for normal operation, please proceed according to the chapter "Initial commissioning".

Malfunctions and damage may only be rectified by authorised trained personnel.

Use only original GREGGERSEN service-parts for repairs!

Service-Parts


No.	Description
325328	Maintenance-kit MC2050„all in one“
325329	Maintenance-kit MC2100„all in one“
325347	Maintenance-kit reserve panel „all in one“
325348	Maintenance-kit Tank reducer panel „all in one“
325258	Maintenance kit for low-pressure reducer MC2050 by date 09. 2001
325260	Maintenance kit for low-pressure reducer MC2050 from date 09. 2001
325257	Maintenance for high-pressure reducer MC2050
325259	Maintenance for high-pressure reducer MC2100
325270	Maintenance kit pneumatic valve assembly
325275	Maintenance kit for connections MC2100
325276	Maintenance kit for connections MC2050
104953	Lithium Battery 3V, 280mAh, CR2430
900859	Programming cable and software Aeolus / Ventus

Electronic switch-over component Aeolus

The electronic switch-over component Aeolus controls and monitors the central gas distribution system. It is located in a separate housing within the system. On the front of the housing there is a lighted display with three LEDs and three keys. The current status of the system is shown on the display. This includes among others the status of each gas source (OPERATION, READY, EMPTY) together with

fault and alarm messages. In the event of an alarm message (an empty message from a supply source is not a fault and therefore does not cause an alarm), in addition a red LED flashes. A window in the hood makes it possible to check the condition of the system on the display even when the hood is closed. Pressure sensors measure the pressure of the cylinder batteries, the tank and the network pressure. Two pneumatic valves (PV1, PV2) control the gas flow of the two cylinder batteries. Potential-free contacts enable the connection of switching electronics to a central monitoring system. A integrated Logbook saves all important system events with time and date. Optional software

enables easier programming of the electronics via the RS485 interface. The most important settings can also be programmed directly at the electronic without a software.

Pos Description

- | | |
|---|---|
| 1 | Network pressure [bar or kPa] |
| 2 | Pressure left side [bar or kPa] (source 1) |
| 3 | Pressure right side [bar or kPa] (source 2) |
| 4 | Tank pressure [bar or kPa] (source 3) |
| 5 | The dashed line show the running source |

When the system is assembled in the factory, the electronic component is programmed to the requirements. The following parameters can be programmed:

- Switch-over between the three versions
 Tank + 2 reserve supply sources (cylinder batteries)
 2 main supply sources + 1 reserve supply
 2 cylinder batteries as reserve supply
- Types of gas
 O₂, N₂O, CO₂, AIR, N₂, SYN, HE, H₂, AR, Misc
- Language of the system messages
 English, German, Polish
- Standard programming
 cylinder batteries: Sensor 0...250 bar, empty 800 kPa, full 2000 kPa
 Tank: Sensor 0...16 bar, empty 700 kPa, full 800 kPa
 Net-pressure: Sensor 0...16 bar, „to low“ 450 kPa, „to high“ 550 kPa

- Inverting the potential-free contacts
- Selection pressure switch/sensor

A mixed use of pressure switches and sensors is not possible. When using pressure switches, the pressure is no longer shown in the display.

If you press the „SET“ button the main menu will be opened.

The program menu is completely in english. In the point of menu “programmed values” the programmed values can be checked. The point

of menu “documentation” shows the saved system events. In the point of menu “system settings” the programmed values (operation mode, gas, configuration source, configuration, outlet, relay control, date/time, language, activation network, reset to factory settings) can be changed.

The current status of the system can be queried via potential-free contacts (relays).

The contacts can be inverted in the settings.

The following potential-free outputs are available:

- left side empty
- right side empty
- tank empty
- network pressure too low
- network pressure too high
- collective Fault
(sensor fault, valve fault, electronic fault, all sources empty, network pressure to high, network pressure to low)
- left side operation
- right side operation
- tank operation

Relay contact is open when the message applies.

Error Message


All alarm messages are shown on the display in plain text. The display changes approx. every 5 seconds from showing the messages to showing the system status. Alarm message (priority 1 and 2) must be confirm from the user. With the “SET” Button you get into the error view and you can choose a single error with the arrow keys. If you press the “SET” Button again the chosen error will be confirm. All error messages will be saved in the documentation (with date and time).

Message display	Priority	LED	Confirm
Sensor 1 defeckt	1	flashing red	to confirm from user (after solve the error)
Sensor 2 defect	1	flashing red	to confirm from user (after solve the error)
Sensor 3 defect	1	flashing red	to confirm from user (after solve the error)
Sensor grid defect	1	flashing red	to confirm from user (after solve the error)
Valve 1 defekt	2	flashing red	to confirm from user (after solve the error)
Valve 2 defekt	2	flashing red	to confirm from user (after solve the error)
Comm error network	3	rot on	automatic after 2s (after fault is solved)
Battery empty	3	rot on	automatic after 2s (after fault is solved)
No gas supply	1	flashing red	to confirm from user (after solve the error)
Netzdruck to high	1	flashing red	to confirm from user (after solve the error)
Network to low	1	flashing red	to confirm from user (after solve the error)
LED defect	3	rot on	automatic after 2s (after fault is solved)

Additional Notes

- An empty message from a supply source is not a fault and therefore does not cause an alarm. It will be virtualized by a flashing yellow LED.
- The battery must be replaced if the electronics indicate it to be necessary. Only authorized technicians are allowed to undertake replacement.

Explosion drawing

Legend explosion drawing

Abbreviation

AN

AVL

AVR

DM1

DM2

DM3

DM4

DS1

DS2

DS3

Assembly

Outgoing to the network

Shut-off valve left

Shut-off valve right

Left high pressure reducer with pneumatic valve

Right high pressure reducer with pneumatic valve

Left low reducer

right low reducer

Pressure sensor left cylinder battery, 0-25000 kPa

Pressure sensor right cylinder battery, 0-25000 kPa

Pressure sensor network pressure, 0-1600 kPa

ELK	Electronic Aeolus
EL	Incoming, left cylinder battery
ER	Incoming, right cylinder battery
HV	Main shut-off valve
M1	Left high pressure gauge
M2	Right high pressure gauge
M3	Medium pressure gauge
M4	Network pressure gauge
NE1	Emergency supply point (NIST)
SVMD	Safety valve medium pressure, 1100 kPa
SVND	Safety valve medium pressure, 700 kPa
ZMD	Connection medium pressure

Plan of terminal connections

The housing of the electronic switch-over component has four connector strips with the following terminals.

Connector 1: 3 pole

Connector 2: 6 pole

Connector 3: 11 pole

Connector 4: 9 pole

Connector 5: 10 pole

Additional Notes:

- Numbering of the plugs from left to right.
- Caution: If you use contact switches you first have to remove the jumpers (W1 ... W4) on the circuit board.
- The contact of the relays is open if the condition is reached.

Connection	Terminal designation	Connector	PI N
Power supply			
24 V DC / AC	+24V / PE / oV	1	1/2/3
Programming interface	RS485- / RS485+ / GND	3	1/2/3
Sensor / contact switch			
Network pressure	+24V / S-Grid	3	4/5
Tank	+24V / S-Q3	3	6/7
HP right side	+24V / S-Q2	3	8/9
HP left side	+24V / S-Q1	3	10/11
Magnetic valves			
Left side	+VQ1 / -VQ1	2	1/2
Right side	+VQ2 / -VQ2	2	3/4
Tank	+VQ3 / -VQ3	2	5/6
potential-free contacts			
Network pressure to low	R1	4	1/2
Network pressure to high	R2	4	3/4
Tank in operation (Q3)	R3	4	5/6
Collective fault	R4-NC / R4-NO / R4-NI	4	7/8/9
Tank empty (Q3)	R5	5	1/2
Right Side empty (Q2)	R6	5	3/4
Left Side empty (Q1)	R7	5	5/6
Right Side in operation (Q2)	R8	5	7/8
Left Side in operation (Q1)	R9	5	9/10

Example electric circuit relay (NI - NC)

Current condition on electronic	Relay contact (NI - NC)
currentless	open
Normal operation	close
alarm	open

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