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

Testing kit for terminal units

Dear customers,
thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0; sales@greggersen.de

Please read the instructions for use thoroughly before using these devices. Only use it once you have completely understood the application and function!

1. Intended use

With this testing kit all necessary measuring- and testing equipment for terminal units for medical compressed gases and vacuum according to EN ISO 7396-1 and -2 is provided. The operation in medical institutions with central gas supply occurs through specialist companies during installation, service and maintenance.

This product is not a medical device. The responsibility for the proper use of the test equipment, the documentation and the interpretation of the measurement results lies with the executing personnel.

2. Product description

The test case consists of the following components:

The test case is equipped with test plugs according to the selected country standard.

1. 902110 Testblock compressed gases for terminal units with flow tube

2. 902115 Testblock vacuum

3. 902116 Vacuum flow tube

4. 326352 Adapter hose

5. 326519 Throttle 25 l with hose

6. 326517 Throttle 50 l with hose

7. 325468 Plug O₂ DIN

- 902130 Plug O₂ BS

- 500219 Plug AIR / O₂ AGA

8. 325466 Plug AIR DIN

- 902131 Plug AIR BS

9. 325469 Plug VAC DIN

- 902132 Plug VAC BS

- 500222 Plug VAC AGA

10. 325470 Plug CO₂ DIN

- 500221 Plug CO₂ AGA

11. 325467 Plug N₂O DIN

- 500220 Plug N₂O AGA

- 902133 Plug N₂O BS

- 902136 Plug AIR 800 BS

- 500268 Plug AIR 800 AGA

- 902136 Plug O₂ / N₂O BS

12. 901261 Leakage testing device DIN


- 902135 Leakage testing device BS

- 902119 Leakage testing device AGA

13. 902112 Particle testing device

14. 103319 Filter paper 00M 55mm

15. 902014 Testing device for Airmotor

Optional components

- 902149 Oxygen-measuring device G 1690-35

- 104478 Substitute-sensor element GOEL370

3. Warnings

- Before each use, please check the components for external damage or contamination.
- Always keep all components free of oil and grease! Danger of fire!
- Pressure fluctuations may occur due to temperature changes. If necessary, these fluctuations in measurements must be corrected according to the equation of state of ideal gases. For this see also attachment E of DIN EN ISO 7396-1.

4. Testing equipment

With the testing kit you have measuring and testing equipment to check terminal units according to DIN EN ISO 7396-1 und 7396-2 and Airmotor-terminal units. These instructions for use explain the individual tests required by normative standards and explain how the measuring and testing equipment works.

The interpretation and documentation of the data obtained with the measuring and testing equipment can be carried out on the basis of DIN EN ISO 7396-1:2019-6.

5. Calibration of the testing equipment

The testing equipment from the test case provided by Greggersen is calibrated by default during production. The resulting calibration protocol is included. Please note that a calibration can only represent a snapshot of the condition of the measuring equipment. It is therefore essential to visually inspect the components and check the plausibility of the measurement results every time they are used. In order to be able to evaluate your tests of systems according to DIN EN ISO 7396-1 and DIN EN ISO 7396-2, in addition to a general visual inspection we recommend every time before usage, to send the test kits to Greggersen Gasetechnik for calibration once a year.

After the calibration has been successfully completed, you will receive a detailed calibration protocol for the individual test equipment from us. In addition, a test sticker with an annual test interval is attached to each test equipment.

6. Tests

The individual tests, as well as their normative origins, that can be carried out with the Greggersen test kit are listed below.

We recommend tests 6.1 to 6.6 for the annual testing of terminal units.

6.1 Leak test

This test is not required by standards, but we recommend it for all terminal units. The Greggersen terminal unit Forano must be carried out before commissioning and during the annual function and leak check.

The following equipment is required for the test:

Leakage testing device

To check the tightness of the terminal unit plug the device into the terminal unit which has to be checked.

This should create a constant dynamic pressure. A leak can be assumed if the pressure drops or increases. If the pressure increases or drops below vacuum, the device must be quickly removed from the terminal unit. The device is not protected against over or under pressure and can be damaged. The leak tester is equipped with a universal plug, making the device suitable for any terminal unit. Possible causes of leaks could be damage to the O-Ring sealings, which then need to be replaced. Furthermore, the closing cone system could be contaminated. At the compressed gas terminal units the pressure on the pressure gauge increases if there is a leak between the closing cone and the lower seal. If the pressure falls, there is a leak in the outer O-ring seal that seals the plug to the terminal unit. At vacuum terminal units the pressure drops if there is a leak in one of the O-Ring sealings. If the pressure drops into the negative, there is a leak in the valve system.

6.2 Test for blockage and flowrate for medical compressed gases and vacuum

Acc. to 12.6.4 / C.3.4.1 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Testblock compr. gases / Vacuum, Vacuum tube, adapter hose, plugs


The test block for compressed gases respectively the test block for vacuum with the vacuum tube must be coupled into the individual terminal units using the plugs. The flow must be adjusted to the value indicated in the table, checking whether the pressure drop exceeds the specified values. If the terminal units are not in a horizontal orientation, the adapter hose must also be used and the measuring tube must be aligned vertically.

Pipe system	Pressure drop	Volume flow
Medical compressed gases other than air or nitrogen for operating surgical tools	$\leq 10\%$	40 l/min
Vacuum	$< 60 \text{ kPa (0,6 bar)}$ (absolute pressure)	25 l/min

6.3 Function- and gas specificity test

This test is not required by standards, but we recommend it for all terminal units. At the Greggersen terminal unit Forano this test must be carried out before commissioning and during the annual function and leak check.

The following equipment is required for the test:

Testblock compr. gases / vacuum, vacuum tube, throttle 25l/min, plugs

The test block for compressed gases respectively the test block for vacuum has to be coupled to the individual terminal unit together with the plugs.

It has to be checked whether:

- The connector is easy to insert, locks and releases properly
- A plug of another type of gas cannot be inserted, locked there or releases gas

6.4 Test for contamination of the pipe manifold system with particles

Acc.to 12.6.10 / C.3.11 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Particle testing device, filter paper, plugs

The housing of the particle tester must be opened and a filter paper must be inserted. Close the particle tester again. The particle testing device must be connected to the terminal unit using the test plug which is suitable for the respective terminal unit so that a gas flow is achieved. The device must be disengaged after 15 seconds. The test device must be opened again and the filter paper must be examined under good lighting conditions. There must be no visible contamination. The filter paper can be reused as long as it is undamaged and as long as no particles from previous tests are visible.

6.5 Test for flow rate and pressure of AGSS-terminal units Type 1L

Acc. to 8.2.2.2 and 12.4.6 of DIN EN ISO 7396-2:2007

The following equipment is required for the test:

Throttle 25l/min and 50l/min, vacuum tube


The terminal unit must be put into operation. The throttles must be connected one after the other to the vacuum flow tube and coupled into the terminal unit. The respective volume flow can be read on the flow tube. With the throttle 50 l/min, the measured flow must not exceed 50 l/min. With the throttle 25 l/min, the measured flow must be at least 25 l/min.

6.6 Test of Airmotor-terminal units for primary- and dynamic pressure

This test is not required by standards, but we recommend it for all Airmotor terminal units. At the Greggersen Airmotor terminal unit this test must be carried out before commissioning and during the annual function and leak check.

The following equipment is required for the test:


Airmotor testing device

During the test, the airmotor terminal unit is checked for the function of the plug connection, the correct pre-pressure and the dynamic pressure in the disposal line. Put the airmotor terminal unit into operation. The test device must be engaged to check the effective connection. The aperture in the test device generates an air flow of 350 l/min. The supply pressure to be read on the upper pressure gauge must be between 700 and 1000 kPa (7 and 10 bar). The dynamic pressure of the exhaust air line shown on the lower pressure gauge must not exceed 20 kPa (0.2 bar).

6.7 Test of area shut off valves for leakage and closure and check for proper division of areas

Acc. to 12.6.2.1 / C.3.2.2.1, C.3.2.2.2, C.3.2.2.3 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Testblock compressed gases, plugs

The test must be carried out after the tests in accordance with C.3.1 have been completed with a positive result. This test is only to be carried out for compressed gases. The system behind an area shut-off valve must be brought to the nominal distribution network. The system must then be isolated from the source at the area shut-off valve. The test block for compressed gases with an adapter plug corresponding the gas type must be connected to the terminal unit. By opening the handwheel, the pressure has to be reduced to 100 kPa (1 bar). The pressure must be noted. The pressure increase should be measured after 15 minutes. A pressure increase of max. 5 kPa (0.05 bar) / 15 min is permitted. If the clarity of the statement is unclear due to the precision of the available measuring instruments, the timeframe can be extended. Consequently, the permissible pressure increase enlarges.

6.8 Test for gas specificity

Acc. to 12.6.16/ C.3.16 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Testblock for compr. gases / vacuum, vacuum tube, oxygen measuring device, plug


All terminal units must be checked using one of the following methods:

- For gases with a specific oxygen concentration, the compressed gas test block must be coupled into the terminal unit together with the respective test plug and a volume flow must be generated. The oxygen sensor must be held in the gas flow. The measurement result must be compared with the expected result (O₂ 100%, O₂ 93%, AIR 21%).
- For gas types without oxygen content, each system must be set to its own pressure value. The compressed gas test block must be coupled into the system together with the respective test plug. The pressure value must be read on the manometer and compared with the set pressure value of the respective gas.
- For vacuum terminal units, the vacuum test block must be coupled into the vacuum terminal unit together with the vacuum test plug. It is necessary to check whether the pressure value displayed on the pressure gauge corresponds to expectations.

6.9 Test for mechanical integrity of vacuum-pipesystems

Acc. to 12.6.1.1 / C.3.1.1 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:


Test block compressed gases, plug VAC

The test must be carried out for the entire vacuum pipe system. The section to be tested must be filled with a test gas with a pressure of 500 kPa (5 bar) and thereafter separated from the source. Precautions must be taken to ensure that the vacuum network is not damaged due to excess pressure. In particular vacuum pressure gauges must be protected from excess pressure or removed before this test. The test block for compressed gases must be coupled to a terminal unit in conjunction with the vacuum plug adapter, whereby the valve must be kept closed. After 5 minutes, check whether there has been a significant drop in pressure in the system and thus a break in the pipe.

6.10 Combined test for leakage and mechanical integrity of vacuum-pipesystems

Acc. to 12.6.1.7 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Test block for compressed gases, plug VAC


The test must be carried out in accordance with 6.2 -> Test for mechanical integrity of vacuum pipe systems.

A possible pressure drop must then be measured over a period of 2 to 24 hours using the same test equipment and the same test pressures. The pressure drop must not exceed 0.025% per hour.

6.11 Test for medical integrity of pipe systems for medical compressed gases

Acc. to 12.6.1.3 / C.3.1.3 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:


Test block compressed gases, plugs

The section to be tested must be filled with a test gas with a pressure of 1.2 times the pressure specified for the section. The test block for compressed gases must be coupled to a terminal unit in conjunction with the plug adapter, whereby the valve must be kept closed. After 5 minutes, check whether there has been a significant drop in pressure in the system and thus a break in the pipe.

6.12 Test for leakage from pipesystems for medical compressed gases

Acc. to 12.6.1.4 / C.3.1.4 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Test block compressed gases, plugs

This test of the compressed gas systems can be carried out in any number of tests of individual sections separated by shut-off valves, whereby no section may be omitted. However, the system must be isolated at least at area shut-off valves and pressure reducers (by setting to zero flow).

The system must be filled with a test gas at the nominal distribution network pressure. The system must then be disconnected from the gas source. The compressed gas test block in combination with the plugs must be coupled into a terminal unit, whereby the valve must be kept closed. Alternatively, the pressure gauge from an area shut-off valve box or a pressure gauge elsewhere in the system can be used.

The pressure and the room temperature must be noted.

After a test period of 2 to 24 hours, the pressure and temperature should be noted again.

The following pressure drops are permitted for this test:


upstream of the area shut-off valve	0,025% / h
downstream of the area shut-off valve	0,4% / h
in areas with flexible hoses	0,6% / h

Note: A change in temperature leads to a pressure change of approximately 0.35%/K and can significantly influence the measurement result of this test. It is recommended to choose a long test period. If in doubt, the test should be repeated at a time when temperatures are more stable.

6.13 Combined test for leakage and mechanic integrity from pipesystems for compressed gases (before covering)

Acc. to 12.6.1.5 / C.3.1.5 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Test block compressed gases, plugs

The section to be tested must be filled with a test gas with a pressure of 1.2 times the pressure specified for the section. The test block for compressed gases must be coupled to a terminal unit in conjunction with the plug adapter, whereby the valve must be kept closed. After 5 minutes, check whether there has been a significant drop in pressure in the system and thus a break in the pipe. The system must be filled with a test gas at the nominal distribution network pressure. The system must then be disconnected from the gas source. The compressed gas test block in combination with the plugs must be coupled into an terminal unit, whereby the valve must be kept closed. Alternatively, the pressure gauge from an area shut-off valve box or a pressure gauge elsewhere in the system can be used.

The pressure and the room temperature must be noted.

After a test period of 2 to 24 hours, the pressure and temperature should be noted again.


A pressure drop of 0.025%/h is permitted for this test.

Note: A change in temperature leads to a pressure change of approximately 0.35%/K and can significantly influence the measurement result of this test. It is recommended to choose a long test period. If in doubt, the test should be repeated at a time when temperatures are more stable.

6.14 Combined test for leakage and mechanic integrity from pipesystems for compressed gases (after covering)

Acc. to 12.6.1.6 / C.3.1.6 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Test block compressed gases, plugs


The section to be tested must be filled with a test gas with a pressure of 1.2 times the pressure specified for the section. The test block for compressed gases must be coupled to a terminal unit in conjunction with the plug adapter, whereby the valve must be kept closed. After 5 minutes, check whether there has been a significant drop in pressure in the system and thus a break in the pipe.

The test for leakage from pipe systems for medical compressed gases must then be carried out in accordance with 12.6.1.4 / C.3.1.4 of DIN EN ISO 7396-1:2019-6.

6.15 Test for leakage in the vacuum pipesystem

Acc. to 12.6.1.2 / C.3.1.2 of DIN EN ISO 7396-1:2019-6

The following equipment is required for the test:

Testblock vacuum, plug VAC

The system should be evacuated until the nominal distribution network pressure is reached. The system must then be separated from the vacuum source.


The vacuum test block, in combination with the vacuum plug adapter, must be coupled into a vacuum terminal unit, whereby the valve must be kept closed. The pressure must be noted. The pressure should be measured again after one hour.

A maximum pressure increase of 20 kPa (0.2 bar) is permitted.

6.16 Test of interconnection

Acc. to 12.6.3 / C.3.3 of DIN EN ISO 7396-1:2019, 12.4.4 of EN ISO 7396-2:2007

The following equipment is required for the test:

Testblock compressed gases / vacuum, vakuum tube, throttle 25l/min, plugs

The sections to be tested must be separated from the supply and the pressure must be relieved. The test block for compressed gases or vacuum must be coupled into the terminal unit together with the respective test plug. Check on the pressure gauges that the entire to be tested area is relieved of pressure. At AGSS terminal units the throttle has to be coupled in. The strand must then be filled with a type of gas. The test equipment must be used to check whether pressure is only present at the terminal unit of this type of gas or whether the AGSS terminal unit is functional and the terminal units lock properly. The test must be carried out at all sampling points for this type of gas. The test must then be carried out for the other types of gas.