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
Ventus 22

Dear customers,

thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggersen.de

Use

For selective removal and monitoring of individual wards / areas from the central gas supply system with integrated pressure monitoring by an integrated emergency alarm unit.

Design


The front is white coated (RAL 9016). The housing of the area valve service unit (AVSU) contains up to 6 area shutoff valves. Each area shutoff valve is equipped with a pressure gauge and a sensor for pressure monitoring. The sensors are internally connected with the area alarm indicator panel. This unit indicates the function of the pipeline system by means of visual indicators. For pressurised gases it is possible to isolate the ward (area) from the supply network physically. For emergency or maintenance reasons, the valve unit is equipped with a NIST- connector. (Optional a second one for HTM guideline).

Operation

In normal operation the door is closed and can be opened with a key. In case of emergency the door can then be unlocked by pressing the key button or push the window inside the valve box. The valve lever direction:

Valve lever vertical in direction of flow = area shutoff valve is open Valve lever horizontal = area shutoff valve is closed.

The display unit together with the integrated area alarm panel is a part of the wall frame and has to clear illuminated panels to indicate normal operational status (green) and alarm conditions for high or low pressure (red), as well as an audible signal for each gas.

Technical data

Design:	Standard in all-metal design, housing of sheet metal, coated white (RAL 9016). One to six area shutoff valves per shutoff valve box with pressure gauge on the valve block (brass), physical separation, emergency supply capability (NIST) and electronic sensor. Inlet and outlet upwards, integrated alarm system.
Operating pressure:	Pressurised gases 400...500 kPa (4...5 bar), max. 1000 kPa (max. 10 bar) Vacuum -99 kPa (-0,99 bar)
Test pressure:	Pressurised gases max. 1500 kPa (15 bar) Vacuum max. 500 kPa (5 bar)
Installation dimensions:	Flush mounted / cavity wall housing Housing 1-3 fold (HxWxD) approx.: 470 x 330 x 77 [mm] Cover 1-3 fold (HxWxD) approx.: 530 x 390 x 16 [mm] Housing 4-6 fold (HxWxD) approx.: 470 x 570 x 77 [mm] Cover 4-6 fold (HxWxD) approx.: 530 x 630 x 16 [mm]
Input / output:	Copper tube dia Ø22x1 [mm]
Alarm system:	Supply voltage: 12...20 V AC oder 12...24V DC
Sensors:	measurement range vacuum / pressurised gas: -100...+600 kPa (-1...6 bar) pressurised gas: 0...1600 kPa (0...16 bar) supply voltage: 10...30 V DC signal: 4...20 mA
Alarm limits: (Standard)	Pressurised gases: $\pm 20\%$ Vacuum: min: -40 kPa; max: -120 kPa
Pressure gauge:	Ø50mm Measurement range: Pressurised gases: 0...1600 kPa (0...16 bar) Vacuum: -100...0 kPa (-1...0 bar)
Guidelines:	DIN EN ISO 7396-1 Medical Devices Directive 93/42/EEC; Classification: IIa

Information

- Floor shutoff valve boxes are part of a central gas supply system according to DIN EN ISO 7396-1. The connecting pipes must be connected by soldering to the gas supply network.
- Comply with the German statutory accident prevention regulations of the BGR 500 Guidelines (UVV Gases / UVV Oxygen)!
- Keep all connections permanently free from oil and grease!
- No warranty claims are accepted in case of damage or malfunction resulting from incorrect transportation, improper storage or incorrect use.
- The product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation.

Recommended installation sequence

- Install the valve units in the back part of the valve box according the gas sequence.
- Install the rear section in / on the wall.
- Connect the shutoff valve box to the supply network (mechanical, electrical).
- After the pressure test, fit the vacuum pressure gauge.
- Connect the sensor cable to the alarm system.
- Connect the cable for further transmission (if present).

! Connect in power-off condition !!!

- Connect the cable of the power supply unit to the alarm system.

! Connect in power-off condition !!!

- Push the front section into the rear section and secure with the 4 screws.
- Open the area shutoff valves and check the line pressure.
- Switch on the power supply for the power supply unit.
- Perform a complete function check according ISO 7396-1 and/or HTM 02-01(2007).

! Check electric security

Installation and start to use is only allowed by certified specialist staff (electric expert). Check the product before the first start-up (initial commissioning) and before every relaunch (after maintenance) according the German Institute for Standardization EN 62353 or Country-specific regulations. Check the permanent connection of the protection conductor (GND).

Measurement of the protection conductor (GND):

1. Measurement: Electronic plate against earth;
2. Measurement: Case against earth;

The limit value is ≤ 0.3 ohms (according to German Institute for Standardization EN 62353, VDE 0751:2008;

- Attach labels to indicate the direction of flow and the type of gas to the left-hand branch of the valve blocks. For compressed gases, the label must be attached so that the arrow points upwards; for vacuum, it must point downwards.
- Perform a complete function check.

Maintenance / service Every shutoff valve box must undergo a regular safety inspection (external damage / pressure values on pressure gauge). Perform an operational test at least once annually on every floor shutoff valve box (function of shutoff valves, pressure gauges / vacuum gauges and alarm system). Malfunctions and damage may only be rectified by authorised trained personnel.

Use only original GREGGERSEN replacement parts for repairs!

Cleaning

Metal front: wipe with a lightly moistened, soft cloth.

PVC window: Wipe carefully, use only light pressure, using soft cloths and commercial glass cleaning agents. Important: never use cleaning agents with abrasive additives or contaminated cleaning water containing particulate matter !!!


Service parts and accessories

Order number	Description
803572	Copper seal 25-18,5 mm
101229	Pressure sensor –1...6 bar with cable
102130	Pressure sensor 0...16 bar with cable
905296	Repa set Mano Ventus 22 O ₂
905297	Repa set Mano Ventus 22 neutral
905298	Repa set Mano Ventus 22 VAC
101079	O-ring, 5 x 1.6 (pressure gauge)
906595	physical separation Ventus 22 10 pcs.
900859	cable and software Ventus

Structure valve unit

Pos.	Description
1	Upper screw connection (inlet), (option NIST connector)
2	Gas shutoff valve
3	Lower screw connection (insert point for physical separation)
4	Fixing hole
5	NIST connection (emergency supply)
6	Pressure sensor
7	Pressure gauge
8	Screw connection outlet
9	Set screw pressure gauge

Replacing pressure gauge

To replace the pressure gauge, carefully unscrew the set screw (9) with hexagon socket head using a suitable hexagon wrench. Hold the pressure gauge during this manoeuvre since it may be under pressure and because of its smooth, unthreaded shaft may jump out due to the excess pressure.

Now pull the pressure gauge off forwards and replace with a new one. Insert the new pressure gauge as far as the stop and then tighten the set screw.

Replacing gas shut-off valves

First seal off the gas conveying line outside the floor shutoff valve box by actuating the last upstream shutoff valve. Then air bleed and depressurise the sealed off part of the piping system. Dismount the pressure gauge and loose the screw holding the brass valve block.

Loosen the upper and lower screw connection and separate from the relevant shutoff valve. Fit a new valve and assemble in the reverse order. Replace the seal of the lower screw connection with a new one. Finally, restore the flow of gas at all closed valves. Check all connections for leaks!

Installation physical separation

To create the physical separation, first place the shutoff valve in the shutoff position by moving the hand lever to the cross position. Then air bleed and depressurise the pipe line downstream of the valve outside the floor shutoff valve box. Then loosen the lower screw connection. lösen. remove the copper seal and replace it by the grommet. Connect the lower screw connection again.

Important: during these operations ensure that the relevant lines are depressurised!!!

On completion of the maintenance / service activities the separation can be removed again. Use a new copper seal.

Inspection

After completing the above assembly, maintenance or repair activities it is essential to check for correctness of installation. In particular, check for leaks and correct operation of electronic components.

After completing the activities, also check whether the flow of each gas is unhindered on opening all previously closed valves and removing the physical separation.

Emergency signal unit

Acknowledge alarms

The active alarm can be muted by pressing the “Reset” button. This means that the buzzer is deactivated for a period of 12 minutes. The displays / relays remain in alarm status. If the cause of the alarm is still present after 12 minutes, the acoustic alarm is reactivated.

If several alarms are activated simultaneously in a multiple warning situation, e.g. Gas 1 pressure too high, Gas 2 pressure too low, the acoustic alarm can be muted for all active alarms by pressing the “Reset” key.

A new, additional alarm immediately activates the buzzer again.

Muting alarms within the alarm – secondary alarm connection via the serial interface

- Press “Reset” on the emergency alarm unit
The alarm is “muted” both on the emergency alarm unit and on the secondary alarm unit
- Press “Reset” on the secondary alarm unit
The alarm is only muted on the secondary alarm unit

Sequence of test alarm

Pressing the “Test” button initiates the following sequence self-test of the EEPROM and program (FLASH) memory of the alarm micro controller:


1. All LEDs off for one second;
2. All red and green LEDs on for one second;
3. Activation of acoustic signaller;
All red LEDs start flashing;
If this is a designer variant, all background lighting is switched on;
4. After 10 sec. the self-test switches itself off.

Any alarm can interrupt the test sequence at any time.

Electrical connections

Pos.	Description
1	Relay contact
2	Sensor connection:1->24V DC (white); 2->Signal (brown)
3	Supply voltage
4	Door sensor
5	Light-emitting diode lighting
6	Interface programming RS485
7	Relays contact main board fault and door alarm
8	Free
9	Connection "M8" board (Option)

current condition emergency alarm unit	relay contact (NI - NC) (dry contact)
currentless	open
normal operating (green)	close
alarm (red)	open

Example electric circuit relay (NI - NC)

Error messages

No.	Description	Prio	LED red	LED green	LED red	Signal
1	normal operating	0	off	on	off	off
2	test alarm	1	flashing	off	flashing	on
3	pressure to high	10	off	off	flashing	on
4	pressure to low	10	flashing	off	off	on
5	door opened	2	off	flashing	off	on
6	pressure switch alarm	10	flashing	off	flashing	on
7	error in module M1	9	flashing	flashing	flashing	off
8	error in module M2/M3	9	off	off	off	on
9	CRC error in EEPROM M1	7	off	flashing	on	on
10	CRC Error in EEPROM M2/M3	7	on	flashing	off	on
11	communication error M1-→M2/M3	8	off	off	off	on
12	communication error M2/M3-→M1	8	flashing	flashing	flashing	off
13	communication error separate alarm unit	3	on	flashing	on	off
14	short circuit at sensor M2/M3	9	on	flashing	flashing	on
15	Sensor cable break on M2/M3	9	flashing	flashing	an	on
16	loss communication main alarm unit -> separate alarm unit	5	on	flashing	on	on
17	Aloss communication with M8	4	blinken	on	on	on
18	error M8	8	on	on	on	off

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