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INSTRUCTIONS FOR USE

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HUMIDIFIER

CE 0482

Deutsch.....	2
English.....	6
Italiano.....	10
Français.....	14
Nederlands.....	18
Polski.....	22
Română.....	26
Magyar.....	30
Hrvatski.....	34
Norsk.....	38

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Humidifier

Dear customers,

Thank you for purchasing this Greggersen product.

If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739357-0, sales@greggersen.de

Read the instructions for use thoroughly before using the medical device. Only use it once you have completely understood the application and function!

Medical device in accordance with:

93/42/EWG, DIN EN ISO 20789

1. Intended use

Refillable humidifier unit for connection to a flow measuring device for humidifying oxygen or compressed air.

2. Functionality

This humidifier unit is a bubble jet humidifier. As such, the humidifier unit is connected to the output of a flow measuring device for oxygen or compressed air. The gas from the dosing unit passes via the 9/16" connection into the bottle container filled with water. The gas bubbles out of the liquid and passes through the outlet barbed connector to the patient. By bubbling through the water, the gas is humidified with vapour or droplets.

3. Variants

The humidifier unit is available as a standard and as an anti-magnetic model.

4. Warnings

- Prior to using the device, the user must be familiar with the product and the instructions for use, and must have understood them.
- In order to prevent possible infections, the cleaning specifications under 10. Reuse must be observed!
- Oxygen is a highly fire-intensifying gas. Always keep all connections free of oil and grease! When handling, avoid the use of hand cream. Do not smoke or expose to open flame near oxygen equipment! Risk of fire!
- Before use, the product must be checked for proper functioning as well as for external damage, obvious leaks and contamination. Non-functional, damaged, leaky or dirty equipment must not be put into operation.
- The humidifier unit must be cleaned regularly according to the specifications in this instruction, especially before every new patient!
- The water level must be checked regularly, at least twice a day.
- Only operate the humidifier with distilled water!
- The humidifier is only to be used with original Greggersen accessories!
- Only use the humidifier in combination with a suitable flow measuring device!
- This humidifier is not suitable for patients whose upper airways have been bypassed with a tracheal tube or tracheostomy. It does not provide sufficient humidification output to avoid irritation of the airways or drying secretion.
- To wash the humidifier, only use cleaning agents intended for clinical use.
- Our products are partly made of brass with a lead content below 1%. Lead is firmly bound in the material as an alloy component, so no exposure is to be expected.
- All serious incidents must be reported to Greggersen Gasetechnik GmbH and to the authority responsible for the place of residence.

5. Symbols

Product contains lead

Manufacturer

Keep free from oil and grease

Consult the instructions for use

Date of manufacture calendar week - year

6. Drawing and parts list

1. Connection nut
2. Device head
3. Outlet nozzle
4. Bubbler
5. Seal
6. Bottle

7. Users

The product is to be used by healthcare professionals.

8. Starting up

Before starting up each time, an external visual inspection and a functional test have to be carried out. Fill the bottle on the humidifier up to the mark with sterile distilled water and screw it onto the device head. Connect the humidifier to a flowmeter or a pressure reducer with flow regulation using the connection nut. The gas flow has to be established for a short time to check functionality by observing the bubbling. Connect the mask or nasal cannula to the outlet barbed connector depending on the type of treatment. Set the flow rate to the required value and start treatment.

9. Taking out of service

To terminate the application, interrupt the flow from the upstream device. The connections can then be released and the bottle emptied. In the event that treatment is briefly interrupted and subsequently resumed with the same patient, the humidifier may remain connected to the upstream device.

10. Reconditioning

The bottle on the humidifier can be reused up to 20 times, provided it is undamaged and uncontaminated. The bottle head has a projected service life of 10 years. Any further re-use, or the use of damaged or contaminated products, is the responsibility of the user.

10.1 Basic principles

All products must be cleaned, disinfected and sterilised before each use. This particularly applies when using for the first time after delivery, as all products are delivered non-sterile. Effective washing and disinfection is an indispensable prerequisite for effective sterilisation.


As part of your responsibility for the sterility of the products, observe that:

- as a matter of principle, only methods that are sufficiently device and product-specific validated are used for washing, disinfection and sterilisation
- the equipment used (disinfectant, steriliser) is regularly maintained and checked
- the validated parameters are adhered to in each cycle

In addition, observe the legal regulations applicable in your country as well as the hygiene regulations of the healthcare institution. This particularly applies to the various specifications with regard to effective prion inactivation.

10.2 Washing and disinfection

If possible, a mechanical method (disinfectant) is to be used for washing and disinfection. A manual method, even using an ultrasonic bath, owing to its significantly lower effectiveness and reproducibility, should only be used if a mechanical method is not available. Pre-treatment must be carried out in either case. The use of a manual washing and disinfection method must be safeguarded by additional product and method-specific validation under the responsibility of the user.

Pre-treatment

Immediately after the application (within 2 h at the latest), coarse impurities must be removed from the products. Use an ultrasonic bath filled with disinfectant solution for this purpose to achieve germ reduction. The disinfectant should be aldehyde-free (otherwise there may be fixation of blood contamination) and suitable for ultrasonic treatment (no foaming), demonstrate tested disinfecting effectiveness (e.g. VAH (German Association for Applied Hygiene) / DGHM (German Society for Hygiene and Microbiology) or FDA approval or CE marking) as well as good washing performance, be suitable for disinfecting instruments made of metals and plastics, and be compatible with the products (see section Material resistance). In addition, use a soft brush or a clean soft cloth to manually remove contamination, which you only use for this purpose, but never metal brushes or steel wool.

Then rinse all openings of the products three times using a single-use syringe (minimum volume 50 ml, if necessary use a suitable tube section as an intermediate piece).

Please note that the disinfectant used during pre-treatment is only for personal protection and is not a substitute for the disinfection step to be carried out later, after washing has been completed.

Mechanical washing/disinfection (disinfector/WD (washer-disinfector))

Care should be taken when selecting a disinfector that

- the disinfector always has tested effectiveness (e.g. DGHM or FDA approval or CE marking according to DIN EN ISO 15883)
- if possible, a tested program is used (in chemical disinfection, risk of disinfectant residues on the products) for thermal disinfection (Ao value > 3000 or min. 5 min. at 90°C for older devices)
- that the program used is suitable for the products and includes sufficient rinsing cycles
- that only sterile or low microbial contamination (max. 10 microbes/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) is used for rinsing
- that the air used for drying is filtered and
- that the disinfector is regularly maintained and checked

Care should be taken when selecting a cleaning agent that

- is generally suitable for cleaning instruments made of metals and plastics
- if, provided thermal disinfection is not used, a suitable disinfectant with tested efficacy (e.g. VAH/DGHM or FDA approval or CE marking) is additionally used, and that this is compatible with the cleaning agent used, and
- that the chemicals used are compatible with the products (see section Material resistance).
- The water used for cleaning must meet the requirements of DIN EN ISO 15883.

The concentrations of the cleaning agent, as specified by the manufacturer and, if necessary, disinfectant, must be strictly adhered to.

Procedure:

1. Disassemble the products as far as possible (bottle, bottle cap, bottle seal).
2. Place the disassembled products in the disinfector. Make sure that the products do not touch each other. Connect all lumina from the products (insufflation rod, humidifier head connected to the hose nozzle) by pushing the rinsing hoses over the rinsing connection of the disinfector.
3. Start the program.
4. Remove the products from the disinfector at the end of the program.
5. Reassemble the products. Instrument oils must not be used.
6. After removing the products, and after additional post-drying in a clean place (if necessary), check and package the products as soon as possible (see section Checking and packaging).

Proof of the basic suitability of the products for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory using the Miele PG 8535 disinfector and the neodisher MediZym cleaning agent 0.5% (v/v). The procedure described above was taken into account here.

10.3 Checking

After washing or washing/disinfecting, check all products for corrosion, damaged surfaces (especially damaged chrome plating of brass parts), chipping and contamination, and discard damaged products (limitation on the number of reuses, see section Reusability). Products that are still soiled must be washed and disinfected again.

10.4 Packaging

Pack the products for sterilization, storage and transport in disposable sterilization packaging (single packaging) that meet the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- suitability for steam sterilisation (temperature resistance up to at least 137°C (279°F) and sufficient steam permeability)
- adequate protection of the products or sterilisation packaging against mechanical damage

10.5 Sterilisation

- Only the sterilisation methods listed below may be used; no other sterilisation methods are permitted. Steam sterilisation
- Fractionated vacuum method or gravitation method (with sufficient product drying)
 - Steam steriliser according to DIN EN 13060 or DIN EN 285
 - Validated according to DIN EN ISO 17665-1 (previously: DIN EN 554 / ANSI AAMI ISO 11134) (valid IQ/OQ (picking) and product-specific performance assessment (PQ))
 - Maximum sterilisation temperature 121°C (250°F; plus tolerance according to DIN EN ISO 17665)
 - Sterilisation time (exposure time at sterilisation temperature) min. 20 min at 121°C (250°F) or min. 5 min at 132°C (270°F)/134°C (273°F).

The use of the less effective gravitation method is only permissible if the fractionated vacuum method is not available; it requires longer exposure times and must be safeguarded by additional product and method-specific validation under the responsibility of the user. Demonstration of the general suitability of the products for effective steam sterilisation has been provided by an independent accredited test laboratory using the Systec V-150 steam steriliser (Systec GmbH Labor-Systemtechnik, Wettenberg) and the fractionated vacuum method. Typical conditions in hospitals and medical practices, as well as the method described above, were taken into account. The flash sterilisation method is generally not permitted. Also, do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation, or plasma sterilisation.

10.6 Storage

After sterilisation, the products must be stored in the sterilisation packaging in a dry and dust-free place.

10.7 Material resistance

When selecting cleaning and disinfecting agents, ensure that the following constituents are not included:

- organic, mineral and oxidising acids (minimum permissible pH 5.5)
- alkalis (maximum permissible pH value 9.5, neutral/enzymatic cleaners recommended)
- organic solvents (e.g. alcohols, ethers, esters, ketones, benzines)
- oxidizing agents (e.g. hydrogen peroxide)
- halogens (chlorine, iodine, bromine), trichlorethylene/perchloroethylene
- aromatic/halogenated hydrocarbons
- ammonia
- oils

Never clean any products with metal brushes or steel wool.

All products must only be exposed to temperatures not higher than 137°C (279°F)!

11. Disposal

There are no special disposal requirements for this device. The device can be disposed of in household waste. If the device is contaminated, take the appropriate precautions.

12. Technical data

Maximum flow rate:	15 l/min
Humidifier system output:	19 mg H ₂ O/l at 1 l/min 13 mg H ₂ O/l at 4 l/min 9 mg H ₂ O/l at 15 l/min
Inlet connector:	9/16"-18 UNF thread
Outlet connector:	Barbed connector according to EN 13544-2
Ambient conditions:	Storage: -30 ... +50°C
	Use: 5 ... +40°C
Dimensions:	90 x 190 x 55 mm (WxHxD)
Weight:	335 g
Classification according to 93/42/EEC:	class IIa non-invasive active therapeutic medical device

13. Accessories

900922 Bottle for sterile H₂O 0.25 l

There is no entitlement to any warranty claim in the event of damage or malfunction arising from improper transport, improper storage or other than the intended use.

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany