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# Compressed gas regulator

## USE

Compressed gas regulators are used to reduce the network pressure of a central gas supply system. The compressed gas regulator is not a medical product and is therefore not suitable for medical applications.

## Variants

The compressed gas regulator is available as a fixed variant (0.5 - 5.5 bar) or as a continuously variable variant (up to 5.5 bar) for oxygen or compressed air.

## Operation

**Compressed gas regulator fixed:** Plug the regulator into the gas outlet in the parking position. Connect accessories to the pressure reducer. Plug the regulator into the gas outlet in the parking position afterwards. After use disengage the regulator from the gas outlet. The regulator can remain in the outlet in the park position.

**Compressed gas regulator adjustable:** Close control valve, plug the regulator into the gas outlet in the operating position. Connect accessories to the regulator and connect to the unit to be supplied. Set the desired working pressure with the control valve. Close the control valve after use. The regulator can remain in the outlet, preferably in park position.

**Compressed gas regulator fixed:** Plug the gas regulator to the parking position of the gas outlet. Connect the accessories onto the outlet of the regulator. Then engage the compressed gas regulator in the service position. After use, uncouple the gas regulator from the gas outlet.

The compressed gas regulator can remain in the parking position of the gas outlet.

**Compressed gas regulator adjustable:** Close the valve, engage the gas regulator in the operating position of the gas outlet. Connect accessories to the gas regulator and connect to the device to be supplied. Set the required working pressure on the valve. After use, close the valve again. The gas regulator can remain in the gas outlet, preferably in the parking position.

## Warnings

- The equipment and all connections must be free from oil, grease and lubricants, otherwise, there is the risk of fire or explosion!
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand. Always unscrew accessories when the pressure reducer is pressureless.
- Cleaning agents and disinfectants must not enter the unit!
- Fluctuations with pre- or counter-pressure, as well as the site conditions may have influence on the flow rate.
- There is no warranty for correct function in case of improper transport or storage or during operating the pressure reducer not in accordance to this instruction of use.
- The product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation.
- Do not use force to close the fine-metering valve.

## Technical data

Design:	plug-in device in chrome-plated all-metal design with pressure gauge 0-6 bar
Type of gas:	oxygen, compressed air
Inlet pressure:	max. 1400 kPa
Back pressure:	fixed or adjustable 0-550 kPa
Connection:	gas type-specific plug insert according to DIN 13260-2, optional international standard
optional:	hose nozzle 7mm
Storage:	-20 ... + 60 °Celsius
Use:	0 ... + 40 °Celsius

## Maintenance/service

External disinfection with a disposable cloth is sufficient for all chrome-plated metal parts. Clean plastic parts with a common disinfectant solution and follow manufacturers instructions. Let equipment dry completely before using it again. Prior to each use, make a visual check (external damage). A function and leak test must be performed at least once a year. Maintenance must be done every five years (replacement of wear parts). Failures and damages may only be remedied by authorized expert personnel. Do repairs only on the pressureless unit. Use only original **GREGGERSEN** spare parts.

## Spare parts / accessories

900441 maintenance kit for compressed gas regulator  
900610 DIN-Connector 7mm, O<sub>2</sub>  
900614 DIN-Connector 7mm, AIR  
900630 DIN-Connector 7mm, CO<sub>2</sub>

[www.greggersen.de](http://www.greggersen.de)

**GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany**