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
INSTRUCTIONS FOR USE

FALKE

Deutsch	2
English	8
Italiano	14
Français	20
Espaniol	26
Nederlands	32
Română	38
Polski	44
Dansk	50
Norsk	56

Druckminderer Falke

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Pressure reducer Falke

Dear customers,

thank you for purchasing this Greggerson product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggerson.de

Please read the instructions for use thoroughly before using this medical device. Only use it once you have completely understood the application and function.

Foreword

The Falke product from Greggerson Gasetechnik GmbH complies with the requirements of (EU) 2017/745 and is classified as an active and non-invasive class IIb device.

Basis UDI: 4031196FALKELU

Intended use

Devices to distribute or supply medical gases, vacuum and/or data to equipment in diagnostic, therapy or surgery.

Device description

The following accessories can be connected to the pressure reducer:

- Flow outlet 9/16-18 UNF: humidifier or nebulizer unit; flow nozzle for oxygen mask or nasal cannula
- Pressure outlet: appliances intended for connection to a central gas supply system

Warnings

- Before using the device, the user must be familiar with the product and the instructions for use and have understood them.
- The appliance may only be used in accordance with its intended purpose.
- Before each use, the product must be checked for functionality and external damage, obvious leaks and soiling or contamination. Non-functional, damaged, leaking, or contaminated devices must not be put into operation.
- Before each use, ensure that the product is correctly aligned.
- When cleaning, ensure that no liquids enter the product. After cleaning, dry the device thoroughly. Do not immerse the device in liquids. The product is not suitable for sterilization.
- Some of our products are made of brass containing more than 1% lead. However, since the lead is securely bound within the alloy, no exposure is expected.
- Only use CE-marked accessories that are suitable for use in combination with this medical device.
- Observe the maintenance and cleaning instructions in this manual.
- Do not close the openings of safety valves.
- No structural modifications may be made to the product.
- Nothing that does not serve the immediate use of the product may be attached to it.
- The product must be handled with care. Avoid unnecessary damage or stress caused by vibrations, impacts, or falls.
- When administering oxygen, please note that an overdose can negatively impact the patient's health. The correct amount to be administered must be determined by medical personnel.
- Oxygen is a highly flammable gas. Always keep all connections free of oils and greases! Avoid using hand cream when handling. No smoking or lighting fires near oxygen equipment! Fire hazard!
- All serious incidents must be reported to Greggersen Gasetechnik GmbH and the local authorities.

Symbols

Product contains lead

Serial number

Consult the instructions
for use

Manufacturer

Date of manufacture,
Calendar week- year

Medical product

Do not use oil

Operation


Before each use, check the product for proper functionality as well as for external damage, obvious leaks, and dirt. Devices that are not functioning properly, damaged, leaking, or dirty must not be put into operation. Before connecting the pressure reducer to the gas cylinder, check the seal in the hand connection fitting. Set the device to "o" on the handwheel (readable through the sight window). Attach the fitting of the pressure reducer to the cylinder valve and tighten the nut by hand. Connect the accessories. Slowly open the cylinder valve. The pressure gauge shows the pressure of the gas cylinder; green area: OK, red area: cylinder empty (cylinder pressure < 5000 kPa)! Set the desired flow rate on the flow adjustment. Check whether the device produces an audible, tactile, or – when using a humidifier – visible gas flow to ensure that the gas supply is functioning properly. Afterwards, the connection to the patient can be made either via the nozzle, the 9/16-18 UNF connection, or alternatively through a humidifier or nebulizer. Caution! The device does not produce flow between two flow stages. After use, set the flowmeter back to "o". If necessary, disconnect the humidifier or nebulizer from the flowmeter. When using humidifiers, ensure that no liquid enters the device. To connect a device to a pressure port, insert the plug into the quick coupling until it audibly and visibly clicks into place, or connect it with the connector and tighten the securing nut. To end the operation, close the cylinder valve and relieve the pressure in the pressure reducer. To do this, briefly open the flow setting until no more gas escapes. If a device is connected at the pressure outlet, open it to relieve pressure.

To disconnect the device, slide the black ring backward to uncouple the plug, or loosen the securing nut.

To deactivate, release the pressure from the reducer and manually disconnect it from the gas cylinder.

1. High pressure connection
2. Pressure gauge
3. Flow outlet
4. Hose connector
5. Pressure connection (coupling; NIST)
6. Flow setting
7. Flow indicator
8. Safety valve

Place of use

The medical staff determines the appropriate use within professional healthcare settings, such as medical facilities. The device is used in hospitals and medical centers across all areas lacking a centralized gas supply. It is intended for reuse on multiple patients.

The devices are available in the following variants and configurations:

- Gas types: O₂, AIR, N₂O, CO₂
- Only with flow outlet
- Only with pressure outlet
- With both - flow and one or two pressure outlets
- Inlet connections are available in various standards: DIN 477, CGA-V1, SS 367615, BS 341-3, NF E29-650, F-627.
- Pressure outlets conform to standards such as DIN 13260-2, NF S 90-116, SS 87 524 30 (AGA), BS 5682, DIN EN ISO 18082 (NIST), M14x1, and G3/8".
- Extended versions are available for large gas cylinders (50 liters), especially when used with humidifiers or nebulizers.

- Devices with various flow rates:

Flow rate	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min
Standard	0,5	1	1,5	2	3	4	6	8	12	15
Standard+	0,5	1	2	4	6	8	10	12	15	30+
Children	0,1	0,3	0,5	0,6	0,8	1	2	3	4	5
Neonatal	0,1	0,15	0,2	0,25	0,3	0,4	0,5	0,6	0,8	1

Cleaning

It is recommended to clean the device after each use through wipe disinfection, following the instructions provided by the disinfectant manufacturer. Use commercially available disinfectants. During cleaning, ensure that no liquids enter the device. The device must be completely dry before reuse. Do not soak the device in disinfectants, as this could cause damage. Always adhere to the hygiene protocols of the facility.

Disposal

There are no special disposal requirements for this device. It can be disposed of with household waste. However, if the device is contaminated, appropriate precautions should be taken.

Maintenance and Service

Before use, visually inspect the device for any external damage. Conduct a functional and leak test annually. Begin by checking the entire device and the O-ring on the cylinder connection for visible damage. Then, operate the device to check for leaks, especially at the safety valve. Verify the functionality of pressure outlets, such as plug-in couplings, by inserting a suitable counterpart. For devices with a flow outlet, set the handwheel to the flow step, ensuring a visible gas flow at the outlet. Maintenance, including the replacement of wear parts, should be performed every five years. All maintenance and repairs must be carried out exclusively by authorized personnel from Greggersen, using original spare parts.

Technical data

Medical device

according to:

DIN EN ISO 10524-1, (EU) 2017/745

Type:

Piston pressure regulator with content pressure gauge,
All-metal design (chrome-plated brass)

inlet pressure(P1):

max. 20000 kPa

Outlet pressure (P2):

450 kPa $\pm$ 50 kPa

Capacity:

max. 100 l/min

Setting accuracy:

$\pm$ 20 % of the set value;
at 1,5 l/min or less $\pm$ 30 % of the set value

Ambient conditions:

Storage:

-20 ... + 60 °Celsius

Useage:

0 ... + 40 °Celsius

Expected service life:

10 years

Spare parts and accessories

900432 Falke DIN repair kit – wear and maintenance

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany