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USER MANUAL

Flowmeter

FLOWMETER

CE 0482

Deutsch	2
English.....	4
Italiano	6
Français	8
Hrvatski	10
Nederlands.....	12
Română.....	14

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
Medical Flowmeter

Intended use

For dosing oxygen, compressed air and other medical gases; in connection with humidifiers or atomizers for vaporizing and atomizing.


Operation

Close the fine-metering valve of the flowmeter before coupling the connector into the gas outlet. Attach rail-mounted units to the unit rail. For plug connection units plug the unit into the gas outlet, for railmounted units plug in the angle plug. The measuring ball will spring upwards briefly as the operating pressure builds up in the flowmeter.

Always connect accessories manually! Sterilise atomizers and humidifiers completely before using for the first time. Screw hose nozzles and humidifiers below the flowmeter. Fill the bottle with sterile water up to the mark and screw on beneath the bottle mount. Connect accessories to the flowmeter.

Set the desired gas quantity with the fine-metering valve. The measuring ball displays the quantity of gas just removed (read-off point: above the ball). Ensure an adequate water level in the atomizer or humidifier during use.

Close the fine-metering valve again after use. The flowmeter may be left in the outlet, ideally in the stand-by position. If required, detach humidifier or atomizer from the flowmeter, and empty and sterilise bottle.

Additional information

- Flowmeter displays correctly only if the measuring tube is vertical.
- Do not close the fine-metering valve using force: the valve seat may be damaged.
- The measuring tube and plummet are coordinated with one another and may only be replaced together.
- Always keep all connections free from oil and grease!
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand. Always unscrew accessoires when the device is pressureless.
- The flowmeter may only be used for the purpose outlined in this instruction of use and by briefed staff.
- Subject to technical change without notice; other models upon request
- There is no warranty for correct function in case of improper transport or storage or during operating the flowmeter not in accordance to this instruction of use.

Technical data

Type:	Pressure-compensated flowmeter for medical gases, display using plummet principle (measuring ball), chrome-plated, all-metal design, plug connection or rail-mounted, optionally available as double flowmeter
Capacity:	15 l/min, infinitely adjustable in the range of the measuring (different flow ranges available)
Operating pressure:	450 kPa $\pm$ 50 kPa
Accuracy of flow:	$\pm$ 10 % of measuring value; for measuring value $\leq$ 5 l/min $\pm$ 0,5 l/min
Inlet:	gas-type specific plug inserts acc. to DIN 13260-2 optional inserts for e.g. SS 8752430, NF 590-116, BS 5682 NIST (DIN EN ISO 18082)
rail units:	9/16-18 UNF (for standard one-way humidifiers)
Outlet:	storage -30 ... + 50 °Celsius
Environment-condition:	0 ... + 40 °Celsius
use:	DIN EN ISO 15002
Guideline:	Medical Products Directive 93/42/ECC
Classification:	Ila


Maintenance/Service

External disinfection with a cloth is sufficient for all chrome-plated brass parts (with customary cleaner). Sterilize humidifier and atomizer bottle at max. 134° C. Let equipment dry completely before using it again.

Prior to each use, make a visual check of flow meter, humidifier and atomizer (external damages). Maintenance must be done every five years (replacement of wear parts). Failures and damages may only be remedied by authorized expert personnel. Use only original **GREGGERSEN** spare parts!

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany