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INSTRUCTIONS FOR USE

NIEDERDRUCK -

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LOW PRESSURE HOSE

CE 0482

Deutsch	2
English	5
Italiano	8
Français	11
Español	14
Nederlands	17
Română	20
Polski	23
Hrvatski	26
Dansk	29
Ελληνικά	32
Norsk	35

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
Low pressure hose

Dear customers,
Thank you for purchasing our
Greggersen product.
If you have any questions or require
information, please contact our sales
and support team.

+49-(0)40 739357-0,
sales@greggersen.de

Please read the instructions for use thoroughly
before using the medical device. Only use it
once you have completely understood the
application and function!

Medical device in accordance with:

DIN EN ISO 5359
(EU) 2017/745

Connections in accordance with:

DIN 13260-2
DIN EN ISO 18082
BS 5682
NF S90-116
SS 875 24 30 (AGA)
SN ENV 737-6 (Carbamed)
UNI 9507
Medap

Symbols:

Consult the instructions
for use

Medical device

Date of manufacture
clender week year

Manufacturer

Keep free from
oil and grease

Product contains lead

Intended use

Supply of medical devices with medical gases or vacuum using a flexible hose line from a central gas supply system or a pressure reducer.

Functionality

The low-pressure hose line establishes a connection between a gas / vacuum source and a medical device or provides gas / vacuum via a gas outlet. In order to prevent a sudden escape of compressed gases when the plug is disengaged, a reverse flow limitation is integrated into our compressed gas hoses. This ensures a controlled release of the gas remaining in the hose after disengaging. Please note that this function is not available for vacuum hoses.

Variants

The low-pressure hose line is available in the following configurations

- gas types O₂, AIR, VAC, CO₂, N₂O, AIR8, O₂/N₂O
- lengths 1.5m, 3m and 5m. Other lengths up to 8 m on request
- neutral colour or with ISO32 marking
- input connections acc. to DIN, BS, NF, SS, SN, UNI, Medap, NIST
- output connections NIST or DIN coupling
- plug handles in angled or straight form
- other versions on request


Warning

- prior to using the low-pressure hose line, the user must be familiar with the product and the instructions for use and has to have understood them.
- the low-pressure hose line must be checked for external damage before each use
- always keep all connections free of oil and grease! risk of fire!
- no smoking or fire near oxygen equipment
- when taking out of service, always disconnect the input connection first
- only connect or release NIST connections in the depressurised state
- never use tools for connecting. All connections can be released by hand
- never connect low-pressure hose lines in series
- low-pressure hose lines cannot be sterilised
- all serious incidents must be reported to Greggersen Gasetechnik GmbH and to the health authority responsible for the place of residence
- our products are partly made of brass with a lead content of more than 0.1%. As an alloy component, lead is firmly bound in the material so that no exposure is to be expected

Drawing and partlist

1. Input connector
2. Gas type identification
3. Hose
4. NIST connection
5. O-ring seal

Operation


Before starting operation each time, an external visual inspection and a functional test have to be carried out. Connect the low-pressure hose line to the device and then connect the input connector to the gas outlet. To take out of service, first disconnect the connection at the gas outlet. When necessary, the connection on the device side can now be released.

Cleaning

Wipe low-pressure hose lines at least once a month with a surface disinfectant based on the active substances aldehydes, quaternary ammonium compounds or alcohol. After cleaning, the product must be completely dried. The hygiene plans of the health facility must be observed.

Servicing

Low-pressure hose lines must be serviced at least every 5 years.

Spare parts

- 102252 O-ring for NIST connector O₂, AIR, N₂O, VAC, AIR/O₂, AIR8
- 102253 O-ring for NIST connector CO₂

Disposal

There are no special disposal requirements for this device. The device can be disposed of in household waste. If the device is contaminated, take the appropriate precautions.

Technical data

Operating pressure :	1400 kPa
Weight:	1,5 m approx. 410 g, 3,0 m approx. 590 g, 5,0 m approx. 810 g
Operating temperature:	-10 °C - 40 °C
Lifespan:	10 years
Classification according to (EU) 2017-745:	non-active, non-invasive, Class IIa
Basis UDI:	4031196NDSCHLAUCHLTG6K

There is no entitlement to any warranty claim in the event of damage or malfunction arising from improper transport, improper storage or other than the intended use.

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany