
GEBRAUCHSANWEISUNG INSTRUCTIONS FOR USE ISTRUZIONI PER L'USO MODE D'EMPLOI

SKUA

CE 0482

Deutsch	2
English	7
Italiano	12
Français	17

Vakuumregler Skua

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Vacuum regulator Skua

Dear customers,

Thank you for purchasing our Greggersen product.
If you have any questions or require informations,
please contact our sales and support team.

+49-(0)40 739357-0, sales@greggersen.de

Please read the instructions for use thoroughly before
using the medical device. Only use it if you have
completely understood the application and function!
All serious incidents must be reported to Greggersen
Gasetechnik GmbH and to the health authority
responsible for the place of residence.

Medical device in accordance with:

93/42/EWG DIN EN ISO 10079-3 DIN EN ISO 13260-2
DIN EN ISO 19054 DIN EN ISO 18082 Carbamed

Symbols

Consult the
instructions for use

manufacturer

Date of manufacture
calender week year

Keep free from
oil and grease

Intended use

The "Skua" vacuum regulator is a device for metering the vacuum
generated from compressed air provided from the central gas supply. It
serves for aspirating fluids and secretion in the medical field.
Suitable for operation in a health facility. The device is operated via a gas
tapping point and covers the performance ranges of high vacuum / high
volume flow.

The central supply system should at least comply with DIN EN ISO 7396-1.
The trained medical personnel determines the application.

Functionality

The vacuum regulator Skua works according to the Venturi principle and generates the vacuum, operated by compressed air.

After connecting the device to the central gas supply via the gas type and standard-specific interface (DIN, BS, NIST, etc.), the desired vacuum can be set by the user using the hand wheel.

Turning the handwheel moves a ram which, via the generated cross-section to the Venturi nozzle, ensures a gas flow and creates the vacuum. The exhaust gas flow is passed through a silencer and released to the environment.

Variants

The Skua vacuum regulator is available in the following versions:


- Power levels -90 and -30 kPa
- Plug devices according to DIN and Carbamed
- Rail devices for DIN rails and Swiss standards

Plug device

Rail device DIN

Rail device Switzerland

Warnings

- The vacuum regulator can only be used by trained medical personnel.
- Not suitable for outdoor operation and/or rescue operation.
- Not suitable for use for drainage in the low vacuum range (e.g. chest drainage).


- For the operation is an secretion overflow bottle with backflow protection required.
- Fluctuations in the upstream pressure and the ambient conditions can affect the flow rate.
- Over-suctioned devices are irreparable and are not allowed to be used be used.

Drawing and part list

1. Hand weheel
2. Manometer
3. Vacuum nozzle
4. Plug or Rail clamp
5. Ovepressure valve
6. Exhaust air sleeve
7. Exhaust air sleeve

Operation

- Before every use, a functional test and external visual inspection has to be carried out to ensure the proper function of the device. Close the vacuum flow, open the hand wheel valve and check whether the device builds up a vacuum. This can be read from the manometer.
- The device must be switched off to connect to the gas supply (fine metering valve closed).
- Do not operate the fine metering valve with force. The regulator unit can be damaged. The valve is closed as soon as a noticeable resistance builds up on closing the fine metering valve.
- Never use tools for connecting the accessories! All connections can be tightened or released by hand. Connecting and releasing the accessories may only be carried out in the depressurised state.
- The exhaust air holes for the driving gas are not allowed to be covered.
- Covering can lead to a pressure increase in the device and therefore a patient risk.
- To avoid malfunctions or faults in operation, the device can only be operated with the outlet directed downward!

Starting up

Plug-in unit:

Couple the vacuum regulator with the connector into the terminal unit. Observe the gas type marking on the label. Ensure that the device is fully engaged.

Rail-mounted unit:

Lock the vacuum regulator in place with a rail clamp on a rail. Ensure that the rail clamp is fully engaged. Then connect the device with the compressed air source via the NIST connection.

Operation

The required vacuum is set via the hand wheel valve and can be read from the manometer. The suction hose is connected at the vacuum nozzle. The manometer can be turned by 360° for improved visibility. The vacuum builds up at different rates depending on the volume of the downstream devices. The direction of the dispelled air can be individually set with the rotatable exhaust air sleeve.

Taking out of service

After using the device, the fine metering valve has to be closed again. As the vacuum is measured at the outlet, the manometer display does not always go to “zero”. Only if the components connected are relieved (opened) does the display go to “zero”. The vacuum regulator can remain in the terminal unit or can be disconnected. Hold the device firmly while disconnecting from the terminal unit!

Cleaning

After each use it is recommended to clean the device with wipe disinfection. Use commercially available disinfectant for this purpose. The device must be fully dried prior to subsequent use. Do not soak the device with disinfectant! The regulator unit can be damaged.

Maintenance and service

- Maintenance has to be performed every 5 years.
- An inspection with a functional test and a leak test has to be performed at least every 2 years.
- In case of unusually heavy-duty use, shorter servicing and inspection intervals may be required.
- Failures and damages may only be remedied by authorized expert personnel.
- Repairs can only be carried out on the depressurised device.
- Use only original GREGGERSEN spare parts!
- There is no entitlement to any warranty claims in the event of damage or malfunction arising from improper transport, improper storage or other than the intended use.

Technical data

Storage:	-30 °C ... +50°C	
Use:	0 ... +40°C	
Gas type:	Compressed air/AIR	
Prim. pressure:	400 ... 500 kPa	
Inlet:	9/16"-18 UNF with nozzle for vacuum-devices	
Skua -30:	Performance class: medium vacuum ajustable: 0 ... -30 kPa ± 10% Volume flow: > 20 l/min	
Skua -90:	Performance class: high vacuum / high volume flow ajustable: 0 ... -85 kPa ± 10% Volume flow > 20 l/min	
Weight:	Plug device:	715 g
	Rail device DIN:	840 g
	Rail device Switzerland:	895 g
Class:	IIa	

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany