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
USER MANUAL

Vakuumregler Spatz

Vacuum Regulator Sparrow

CE 0482

Deutsch	2
English	5
Italiano	8
Français	11
Nederlands	14
ελληνικά.....	17
Hrvatski	20

Vakuumregler Spatz

Verwendung

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
Vacuum regulator Sparrow

Use

To adjust the vacuum power for absorb medical liquids. The following accessories can be connected to the vacuum regulator:

Accessories

- overflow bottle
- tube-connection 9/16

Operation

Plug-in units can be directly plugged into an outlet; rail units are fastened on the medical instruments rail; vacuum supply is via the hose connection.

Prior to plugging the unit into the outlet, close the needle valve and the quick-action valve of the ejector regulator. For plug-in units, plug the unit into the outlet; and for rail units, plug the angle plug into the outlet.

The vacuum regulator does **not** have any overflow safeguarding!

Mandatory necessary for the mode: Operating equipment always use with a overflow bottle. With secret contaminated equipment is not repairable

Open the quick-action valve completely and adjust the desired vacuum with the needle valve. Suction can be interrupted any time with the quick-action valve without changing the set value.

Turn the valve to the left to adjust the desired vacuum. The safety valve will open when the adjusted vacuum is more then -16 kPa "Sparrow -16" or -10 kPa "Sparrow -10"; a pulsating sound will appear.

After use, close the needle valve and the quick-action valve again. The vacuum regulator can remain in the outlet, preferably in park position.

Additional Information

- The vacuum regulator should only be operated by authorized personnel!
- Spatz -90: the bypass bore must always be free since the fine vacuum adjustment would otherwise be impaire.
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand.
- The vacuum regulator Sparrow -16 and -10 are limited by a safety valve; keep the holes free at all times; modification of the safety valve may influence the safety of the device!
- Do not forcibly close needle valve and quick-action valve: the valve seat may be damaged!
- Always operate the vacuum regulator with overflow safety: liquids or solids in the unit impair its performance!
- Subject to technical change without notice; other designs upon request.

Additional information

- Fluctuations with pre or counter-pressure, as well as the site conditions may have influence on the flow rate.
- The product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation.
- There is no warranty for correct function in case of improper transport or storage or during operating the vacuum regulator not in accordance to this instruction of use.

Technical data

Type:	chromium plated all-metal construction, plug- or rail unit, vacuum adjustment by valve, vacuum gauge, Sparrow -16 and -10 with safety valve	
Capacity: (Sparrow -90)	nominal range pressure gauge: infinitely adjustable between	0 ... -100 kPa -20 ... -85 kPa
Capacity: (High-Sparrow -90)	nominal range pressure gauge: infinitely adjustable between	0 ... -100 kPa 0 ... -100 kPa
Capacity: (Sparrow -16)	nominal range pressure gauge: infinitely adjustable between vacuum is limited by safety valve	0 ... -25 kPa 0 ... -16 kPa (-160 mbar)
Capacity: (Sparrow -10)	nominal range pressure gauge: infinitely adjustable between vacuum is limited by safety valve	0 ... -16 kPa 0 ... -10 kPa (-100 mbar)
Gas:	vacuum	
Operating pressure:	-99 kPa ... -40 kPa	
Inlet:	gas-type specific plug inserts according to DIN 13260-2 optional inserts for e.g. SS 8752430, NF S90-116, BS 5682 rail units: NIST (DIN EN ISO 18082)	
Outlet:	9/16-18 UNF with vacuum tube-connection	
storage:	-30 ... + 50 °Celsius	
use:	0 ... + 40 °Celsius	

Guideline: DIN EN ISO 10079-3

Medical Devices Directive 93/42/ECC; Classification: IIa

Maintenance/Service


External disinfection with a disposable cloth is sufficient for all chrome-plated metal parts. Clean plastic parts with a common disinfectant solution and follow manufacturers instructions. Let equipment dry completely before using it again.

Prior to each use, make a visual check (external damage). A function and leak test must be performed at least once a year. Maintenance must be done every five years (replacement of wear parts). Failures and damages may only be remedied by authorized expert personnel. Do repairs only on the pressureless unit. Use only original **GREGGERSEN** spare parts. The O-Ring must be replaced regularly.

Spare parts /accessories

item no.	Description
900628	tube-connection 9/16" vacuum
900913	secretion overflow bottle, 0,25L

pos	Description
1	vacuum gauge
2	needle valve
3	quick-action valve
4	vacuum nipple with nut
5	safety valve (only Sparrow -16 und -10)

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany