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
INSTRUCTIONS FOR USE

SPECHT

VAKUUMREGLER
VACUUM REGULATOR

CE 0482

Deutsch	2
English.....	5
Italiano	8
Français	11
Nederlands	14
Hrvatski.....	17
Polski.....	20
Dansk	23

Vakuumregler Specht

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
Vacuum regulator Specht

Intended use

To adjust the vacuum power for absorb medical liquids.

The following components can be connected to the vacuum regulator:

- overflow bottle
- tube-connection 6mm

Operation

Plug-in units can be directly plugged into an outlet; rail units are fastened on the medical instruments rail; propellant gas is supplied via the hose connection.

Prior to plugging the unit into the outlet, close the needle valve and the quick-action valve of the vacuum regulator. For plug-in units, plug the unit into the outlet; and for rail units, plug the angle plug into the outlet.

The vacuum regulator does **not** have any overflow safeguarding!

Mandatory necessary for the mode: Operating equipment always use with a overflow bottle. With secret contaminated equipment is not repairable

Open the quick-action valve completely and adjust the desired vacuum with the needle valve. Suction can be interrupted any time with the quick-action valve without changing the set value.

Turn the valve to the left to adjust the desired vacuum. The safety valve will open when the adjusted vacuum is more then -16 kPa "Specht -16" or -10 kPa "Specht -10"; a pulsating sound will appear.

After use, close the needle valve and the quick-action valve again. The vacuum regulator can remain in the outlet, preferably in park position.

Warnings and Notices

- The vacuum regulator should only be operated by authorized personnel!
- Propellant gas must be able to freely escape through the outlet bore since the fine vacuum adjustment would otherwise be impaired. A blocked outlet bore can be result into a positive pressure in the regulator and this could cause to a compromise to the patient.
- The ejector sleeve is conically designed to accept a bacterial filter.
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand.
- The vacuum regulator Specht -16 and -10 are limited by a safety valve; keep the holes free at all times; modification of the safety valve may influence the safety of the device!
- Do not forcibly close needle valve and quick-action valve: the valve seat may be damaged!
- Always operate the vacuum regulator with overflow safety: liquids or solids in the unit impair its performance!
- Subject to technical change without notice; other designs upon request.

Warnings and Notices

- After establishing the connection to the supply network, the device must be put into operation for a short time. This checks that the connection has actually been established and that, for example, the device plug is not in the parking position of the gas outlet.
- Fluctuations with pre or counter-pressure, as well as the site conditions may have influence on the flow rate.
- The product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation.
- There is no warranty for correct function in case of improper transport or storage or during operating the vacuum regulator not in accordance to this instruction of use.
- Regulator must not be used for drainage in the low-vacuum range (e.g. thoracic drainage), central power supply (by generating vacuum/compressed air), pipe systems of vehicles and buildings as well as wall connections, end pieces such as suction catheters, Yankauer suction cups and suction nozzles, syringes, dental suction devices, anesthetic gas Traction Systems, Laboratory Suction, Autotransfusion Systems, Closed Wound Drainage Systems, Mouth Operated Mucus Aspirators including Neonatal Mucus Aspirators, Obstetric Devices, Breast Pumps, Liposuction, Uterine Suction and Smoke Evacuation Systems.
- When starting up, make sure that all screwable connections are tightened hand-tight.
- All variants can be used for pharyngeal suction.
- Only devices marked as „antimagnetic“ are suitable for use in the MRI area.
- The device is not approved for outdoor use and/or transport.
- If liquids or solids have entered the device, it must be disposed of.

Technical data

Type:	chromium plated all-metal construction, plug- or rail unit, vacuum adjustment by valve, vacuum gauge, Specht -16 and -10 with safety valve		
Capacity: (Specht -90)	nominal range pressure gauge:	0 ... -100 kPa	
	infinitely adjustable between	0 ... -90 kPa	
	maximum flow:	> 20 l/min	
Capacity: (Specht -60)	nominal range pressure gauge:	0 ... -100 kPa	
	infinitely adjustable between	0 ... -60 kPa	
	maximum flow:	> 25 l/min	
Capacity: (Specht -16)	nominal range pressure gauge:	0 ... -25 kPa	
	infinitely adjustable between	0 ... -16 kPa (-160 mbar)	
	maximum flow:	> 18 l/min	
Capacity: (Specht -10)	nominal range pressure gauge:	0 ... -16 kPa	
	infinitely adjustable between	0 ... -10 kPa (-100 mbar)	
	maximum flow:	> 18 l/min	
Gas:	compressed air		
Operating pressure:	450 kPa ± 50 kPa		
Inlet:	gas-type specific plug inserts according to DIN 13260-2 optional inserts for e.g. SS 8752430, NF S90-116, BS 5682 rail units: NIST (DIN EN ISO 18082)		
Outlet:	9/16-18 UNF with vacuum tube-connection		
Environment-condition:	storage:	-30 ... + 50 °Celsius	
	use:	0 ... + 40 °Celsius	
Guideline:	DIN EN ISO 10079-3 Medical Devices Directive 93/42/ECC; Classification: IIa		

Maintenance/Service


External disinfection with a disposable cloth is sufficient for all chrome-plated metal parts. Clean plastic parts with a common disinfectant solution and follow manufacturers instructions. Let equipment dry completely before using it again.

Prior to each use, make a visual check (external damage). A function and leak test must be performed at least once a year. Maintenance must be done every five years (replacement of wear parts). Failures and damages may only be remedied by authorized expert personnel. Do repairs only on the pressureless unit. Use only original **GREGGERSEN** spare parts. The O-Ring must be replaced regularly.

Spare parts /components necessary for operation

Artikelnr.	Bezeichnung
900628	tube-connection 9/16" vacuum
900911	suction bottle 0,25l - AIR

pos	description
1	vacuum gauge
2	needle valve
3	quick-action valve
4	filter take-up
5	vacuum nipple with nut
6	safety valve (only Specht -16 und -10)

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-31 / 21031 Hamburg / Germany