Varioport Specht Varioport Woodpecker

Verwendung

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
Use

To adjust the vacuum power for the absorption of medical liquids. The following accessories can be connected to the vacuum regulator:

Accessories

- secretion- and rinse bottle
- tube-connection 9/16
- sleeve

Technical data

Type: chromium plated all-metal construction, rail unit, colour coated, vacuum adjustment by valve, vacuum gauge

Output:	nominal range pressure gauge:	0 ... -100 kPa
(Woodpecker -90)	infinitely adjustable between:	0 ... -90 kPa
	maximum flow:	> 20 l/min

Gas: compressed air

Fuel gas: at -85 kPa $\approx$ 30 l/min, at -40 kPa $\approx$ 15 l/min

Operating pressure: 450 kPa $\pm$ 50 kPa


Inlet: NIST (DIN EN ISO 5359)

Outlet: 9/16-18 UNF with vacuum tube-connection

Environment-condition:	storage:	-30 ... + 50 °Celsius
	use:	0 ... + 40 °Celsius

Guideline: DIN EN ISO 10079-3
DIN EN ISO 5359
DIN EN ISO 19054
Medical Products Directive 93/42/ECC; Classification: IIa

Dimensions for bottles

Additional information

- The vacuum regulator should only be operated by authorized personnel!
- Propellant gas must be able to freely escape through the outlet bore, otherwise the fine vacuum adjustment would be impaired.
- The ejector sleeve is conically designed to accept a bacterial filter.
- Do not forcibly close needle valve and quick-action valve: the valve seat may be damaged!
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand.
- Always operate the vacuum regulator with overflow safety: liquids or solids in the unit influence the performance!
- Fluctuations with pre or counter-pressure, as well as the site conditions may have influence on the flow rate.
- There is no warranty for correct function in case of improper transport or storage or operating the vacuum regulator not in accordance to this instruction of use.
- Subject to technical change without notice; other designs upon request.

Operation

Plug-in units can be directly plugged into an outlet; rail units are fastened on the medical instruments rail; propellant gas is supplied via the hose connection.

Prior to plugging the unit into the outlet, close the needle valve and the quick-action valve of the ejector regulator. For plug-in units, plug the unit into the outlet; and for rail units, plug the angle plug into the outlet.

The vacuum Regulator does **not** have any overflow safeguarding!

Mandatorily necessary for the mode: Always use operating equipment with an overflow bottle. Secret-contaminated equipment is not repairable

Open the quick-action valve completely (turn by 180°) and adjust the required vacuum with the needle valve. Suction can be switched off at any time with the quick-action valve without changing the set value.

After use, close the needle valve and the quick-action valve again.

Maintenance/Service


External disinfection with a cloth is sufficient for all chrome-plated brass parts (with customary cleaner). Equipment needs to dry completely before use.

Prior to each use, make a visual check of regulator (external damages). A function and leak test should be performed at least once a year. The units should have a maintenance check at least once every five years. If the regulator is used to an unusual high degree, the intervals between maintenance might be shorter. Failures and damages may only be repaired by authorized personnel. Use only original **GREGGERSEN** spare parts!

Spare parts / Accessories

number	name
900628	tube-connection 9/16
110048	sterile-filter air-ejector
802694	sleeve Varioport (2 piece)

Pos	Bezeichnung
1	vacuum gauge
2	needle valve
3	quick-action valve
4	filter take-up
5	NIST-connection
6	rail clamp
7	vacuum nipple with nut

Greggersen Gasetechnik GmbH

Bodestrasse 27-29

21031 Hamburg

Germany

Tel.: +49 (0)40 - 73 93 57 - 0

Fax: +49 (0)40 - 73 93 57 - 39

E-Mail: info@Greggersen.de

Internet: www.Greggersen.de