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
# Ventus evo 2


## Dear customers,

Thank you for purchasing this Greggersen product.

If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739357-0, [sales@greggersen.de](mailto:sales@greggersen.de)

Please read the instructions for use thoroughly before using the medical device. Only use it once you have fully understood the application and function!

## 1. General

The Ventus evo 2 product from Greggersen Gasetechnik GmbH complies with the requirements of (EU) 2017/745 and is classified as a Class IIa active, non-invasive device.

### 1.1 Medical device in accordance with:

DIN EN ISO 7396-1

HTM 02-01 Part A

### Configuration group

The Ventus evo 2 is a configurable medical device. The complete zone control and shut-off unit is assembled from several individual components in accordance with the respective requirements. The configuration group comprises all permitted variants of the Ventus evo 2 system. A zone control and shut-off unit consists of a rear part and a front part containing electronics and valve inserts with sensors.

A fully configured and correctly assembled Ventus evo 2 zone control and shut-off unit is authorised as a medical device under the number 902249 Ventus evo 2 configuration group. The UDI marking for the medical device is located on the rear part of the zone control and shut-off unit.

In accordance with HTM 02-01 Part A, valve inserts with a NIST connection are required upstream of the shut-off valve. This requirement is met exclusively by configurations in which valve inserts with two NIST connections are used.

### 1.2 Intended purpose

Devices to distribute or supply medical gases or vacuum to equipment in diagnostic, therapy or surgery

### 1.3 Intended users

Intended users are medically or technically trained professionals.

#### 1.4 Target patient group

All patient groups requiring a supply of medical gases or vacuum, or treatment with gas-powered therapeutic devices. The product is intended to supply gases or vacuum. The decision on its use rests with the medical staff.

#### 1.5 Warnings

- Before using the device, the user must be familiar with the product and the instructions for use, and must have understood their contents.
- The device must only be operated by trained personnel who are familiar with the procedures of the respective healthcare facility.
- Installation, commissioning, maintenance, repair and decommissioning must be carried out exclusively by qualified specialist personnel.
- Prior to commissioning following a new installation, modifications, or maintenance or repair work involving the dismantling of essential gas-carrying components, a commissioning test in accordance with DIN EN ISO 7396-1 must be carried out.
- During installation, the pipework must be purged with inert gas whilst soldering. The housing and adjacent components must be protected from thermal stress.
- Electrical connections and maintenance work on the product, including on gas-carrying components, must only be carried out when the system is de-energised.
- To ensure safe decommissioning, the device must be connected via a dedicated fuse or a separate electrical circuit.
- Zone control and shut-off units must be positioned so that they are clearly visible and freely accessible at all times.
- During installation, cleaning, maintenance and repair, care must be taken to ensure that no contaminants enter the gas-carrying parts.
- Prevent liquids, in particular cleaning agents or leak detection agents, from entering the interior of the unit or coming into contact with electronic components.
- Only original spare parts from Greggersen Gasetechnik GmbH may be used for maintenance and repair work.
- During assembly, the front part, which has not yet been screwed into place, may tip out of the rear part when the door is open. It must therefore be secured until it is fully bolted in place.
- HF communication devices must not be operated or used at a distance of less than 30 cm from the unit.
- All serious incidents relating to the product must be reported immediately to Greggersen Gasetechnik GmbH and the relevant national health authority.



## 2. Unit description

The Ventus evo 2 zone control and shut-off unit is designed for the targeted shut-off and monitoring of individual zones or stations within a medical gas supply system. It enables the safe isolation of a defined supply zone from the central gas network, as well as the continuous monitoring of the gas pressure of the connected media.

The zone control and shut-off unit consists of a powder-coated metal housing with a side-opening door and a stainless steel rear part. The door is available with or without a viewing window and with either a lock or a push-button. In versions with a lock, the door can be opened in an emergency by pressing in the lock cylinder.

1. Valve insert
2. Door
3. Lock / Push-button
4. Window
5. Display
6. Inlet pipe
7. Copper seal
8. Ball valve
9. Outlet pipe
10. Pressure gauge
11. Mounting hole
12. Pressure gauge mounting screw
13. Sensor
14. NIST connection
15. NIST cap

Depending on the model, up to six zone shut-off valves for different medical gases can be integrated. The valve inserts route the respective gas from the source to the supply zone. An integrated ball valve allows the relevant zone to be manually disconnected from the central gas supply system.

Each zone shut-off valve is fitted with a pressure sensor and a pressure gauge. The pressure sensor continuously monitors the gas pressure in the respective zone and transmits the measured values to the integrated zone monitor or the emergency signalling device. There, the readings are compared with the individually set limit values for the healthcare facility. The operating status of the individual media is indicated on the display. In addition, an audible warning is triggered in the event of an alarm. The audible alarm can be acknowledged by authorised personnel.

Optionally, the zone control and shut-off unit can be fitted with a door monitoring system. If the door is opened without authorisation, a visual and audible alarm is triggered. This can be deactivated by trained staff. Before opening the door, the door-opening alarm can be deactivated using a key combination of the test and reset buttons.

The zone control and shut-off unit is fitted as standard with a NIST emergency inlet point downstream of the ball valve. Optionally, an additional NIST emergency inlet point can be

integrated upstream of the ball valve.

These connections enable an external gas supply to the affected area in the event of an emergency.

To forward alarm and status messages to a building management system, the zone monitor can be expanded with up to five relay boards. A potential-free contact is also available. A power failure is indicated by a signal light on the front panel.

The zone control and shut-off unit is available in various sizes for 1 to 3 or 4 to 6 media, and optionally with a surface-mounted frame.

### 3. Installation and Commissioning

Installation and commissioning must only be carried out by qualified specialist personnel. Prior to initial commissioning and any subsequent recommissioning, the product must be inspected in accordance with the applicable national regulations and standards for electrical equipment.

#### 3.1 Checking the delivery

Before installation, the product must be checked for damage to the transport packaging, for obvious damage to the product itself, and to ensure that the delivery is complete and the configuration is correct in accordance with the delivery note.

#### 3.2 Required tools and equipment

To install the Ventus evo 2 zone control and shut-off unit, only standard tools such as a drill, cordless screwdriver, Allen key, spirit level, wire strippers and multimeter are required. A putty knife and filler may be required to seal any gaps.


#### 3.3 Installation

Remove the housing from the packaging and dismantle the front part. To do this, open the door, loosen the screws above the display, fold the electronics forward and remove them by lifting them slightly. Cut out openings in the rear part for electrical cables and for the copper pipes.

##### 3.3.1 Flush-mounted installation in concrete or masonry walls


Chisel out the installation opening for the housing and mount the housing horizontally in the wall using the rear or side holes. The lower front edge of the housing must be flush with the finished wall. The thickness of the plaster layer, tiled panels, etc. must be taken into account here. Fit the valve assemblies in the rear part according to the configuration. Feed the electrical cables into the housing.

Pipe culverts are cut to size during assembly up to this stage.



### 3.3.2 Installation in plasterboard walls

Cut an installation opening in the plasterboard wall. The opening must not exceed the dimensions of the housing by more than 5 mm on any side. Where applicable, take existing fire protection enclosures into account. Secure the housing horizontally to the stud frame of the plasterboard wall using the provided fixing holes. Alternatively, the rear part can be screwed directly to the plasterboard wall using the Ventus evo 2 mounting brackets for drywall installation. The mounting system is designed so that there is a gap of 12.5 mm between the front edge of the mounting system and the top edge of the bracket. Fit the valve assemblies into the rear part according to the configuration. Feed the electrical cables into the housing.

### 3.3.3 Surface-mounted installation

Adhere the spacers to the rear of the rear part box. Then screw the rear part box to the wall through the rear fixing holes. Next, position the cover plates flush against the wall, slide them up against the box and secure them via the slotted holes.

### 3.4 Gas connection

Open the ball valve during the brazing process and purge the pipework with inert gas. Take care to avoid heat damage and, if necessary, protect the casing and pipework. Afterwards, cover the rear part with the protective cover. Before closing up the walls, carry out a pressure test on the system in accordance with DIN EN ISO 7396-1. In the case of vacuum valve inserts, note that the pressure test must be carried out without a pressure gauge, as this would otherwise be damaged. For this purpose, the valve block is fitted with a blind plug on delivery.

Following the pressure test, remove the blind plug and fit the pressure gauge. See 6.2 'Replacing the pressure gauge' for details.

### Applying the labels indicating the direction of flow

After soldering the pipes, the labels indicating the direction of flow and the type of gas must be affixed to the left-hand branch of the valve blocks. For pressurised gases, the label must be applied so that the arrow points upwards; for vacuum, it must point downwards.

### 3.5 Fitting the front part

To fit the front part, slide it onto the rear part and secure it with the screws from the inside. To ensure proper earthing and thus electrical safety, make sure that the toothed washer pre-fitted in the rear part is visible beneath the screw head after assembly.

### 3.6 Electrical connection

The electrical connection must only be carried out when the power is off!

Strip 8–9 mm of insulation from the cable ends. Use ferrule sleeves for flexible conductors.

Recommendation for cable selection

Mains voltage terminal blocks:	NYM – 3x1.5 mm <sup>2</sup> to 3x2.5 mm <sup>2</sup>
Analogue/digital inputs:	LiYCY – 2x2x0.5 mm <sup>2</sup> (flexible) Do not connect the shield IStY – 2x2x0.6 mm <sup>2</sup> (rigid)
Relay modules:	IStY – 2x2x0.6 mm <sup>2</sup> , IStY – 6x2x0.6 mm <sup>2</sup>
Network cable:	CAT 6 or higher

#### 3.6.1 Installation of electrical components

The electrical components of the Ventus evo 2 are supplied pre-assembled as a module. On site, the only task required is to mount this module onto the designated top-hat rails.

The assembly is positioned on the top-hat rails as follows:


- Secure the terminal blocks for the mains connection to the right-hand top-hat rail.
- Mount the power supply unit on the left-hand top-hat rail.
- Screw the mains filter into place centrally.

Once mounted on the top-hat rails, the following connections must be made:

- Connect the 24 V supply cable to the electronics.
- Connect the spare earth cable to the designated earth point on the electronics.

Subsequently, check that all components are correctly secured and that the connection cables are routed safely.

**Note:** Plugging in the mains plug and installing the jumper activates the visual alarm via an LED in the event of a mains power failure. To prevent the internal battery from discharging too quickly, the mains plug must be disconnected or the jumper removed during transport, commissioning or other prolonged periods without a power supply.

#### 3.6.2 Connecting sensors and the relay module

1. Connect the sensor cables to the electronics.

The pin assignment is shown on the rear of the electronics.

2. Secure the mains filter with screws.
3. Clamp the power supply unit onto the DIN rail.
4. Fit the connecting cables.


### 3.6.3 Connecting the base board

Digital 1-4	
Pin	Function
1	Channel 4
2	Channel 4
3	Channel 3
4	Channel 3
5	Channel 2
6	Channel 2
7	Channel 1
8	Channel 1

Digital 5-8	
Pin	Function
1	Channel 8
2	Channel 8
3	Channel 7
4	Channel 7
5	Channel 6
6	Channel 6
7	Channel 5
8	Channel 5

Analogue 1-3	
Pin	Function
1	Channel 3-
2	Channel 3+
3	Channel 2-
4	Channel 2+
5	Channel 1-
6	Channel 1+

Analogue 3-6	
Pin	Function
1	Channel 6-
2	Channel 6+
3	Channel 5-
4	Channel 5+
5	Channel 4-
6	Channel 4+

Supply voltage		
Pin	Function	Note
1	+24 V DC	
2	GND	
3	Battery jumper	See 3.6.1
4		


Relay module		
Pin	Function	Cable colour
1	+24 V DC	red
2	D -	white
3	D +	green
4	GND	black

Relay	
Pin	Function
1	NC
2	C
3	NO


### 3.6.4 Connecting the relay module

#### Connecting the relay module to the base board


#### Addressing the relay modules

RM	S1	S2	S3	S4
1	OFF	OFF	OFF	OFF
2	ON	OFF	OFF	OFF
3	OFF	ON	OFF	OFF
4	ON	ON	OFF	OFF
5	OFF	OFF	ON	OFF

#### Numbering of the potential-free contacts


Relay 1 - 3	
Pin	Function
1	Relay 3
2	Relay 3
3	Relay 2
4	Relay 2
5	Relay 1
6	Relay 1

Relay 4 - 6	
Pin	Function
1	Relay 4
2	Relay 4
3	Relay 5
4	Relay 5
5	Relay 6
6	Relay 6

### 3.7 Plaster protection

To protect the rear part, including any valve inserts that may already be fitted, during the construction phase, Ventus evo 2 products are supplied with a plaster protection cover. This protects the components from dirt and damage caused by plastering, painting and fitting work. The plaster protection must be folded in accordance with the printed instructions and fitted completely over the rear part. It remains on the unit until the finishing works are complete and must only be removed prior to final assembly. This prevents dirt from entering the product.


### 3.8 Optional conversion kit for plug-in pressure sensor

The optional conversion kit for plug-in pressure sensors enables the pressure sensors, the monitoring electronics and the signal transmission to be checked whilst the unit is in operation. To do this, the pressure sensor can be removed from the valve insert by loosening a threaded pin and connected to a suitable test device. An integrated vent port minimises any potential gas leakage during testing.

#### Fitting the conversion kit

1. Unscrew the pressure sensor from the valve insert.
2. Remove the existing seal and dispose of it correctly.
3. Screw the connector fitting, fitted with a new seal, into the valve insert.
4. Position the threaded pin in the connector fitting so that it is accessible from the front.
5. Screw the pressure sensor, fitted with a new seal, into the elbow.
6. Insert the sensor and elbow into the connector and secure with the threaded pin.
7. Once the system has been pressurised, check all connections in the conversion kit for leaks using a suitable leak detection fluid.

## Removing and reinserting the sensor for testing

1. Loosen the set screw and pull the sensor out of the connector.
2. Insert the blind plug into the valve block and secure it with the set screw.
3. Test the sensor in a suitable test rig.
4. Loosen the set screw and remove the blind plug.
5. Reinsert the sensor into the valve block and secure it with the set screw.
6. Check the connections of the conversion kit using leak detection fluid.

## 4. Operation

The Ventus evo 2 zone control and shut-off unit is designed to be intuitive to operate and can be used by untrained personnel in an emergency.

### 4.1 Opening the door

Depending on the model, the door of the Ventus evo 2 can be unlocked and then opened either by operating the lock with the corresponding key or by pressing the push-button. In an emergency – for example, in the event of a gas leak or fire – as well as for maintenance work, it may be necessary to open the door to gain access to the interior of the zone control and shut-off unit. For versions with a push-button, the door can be opened directly for this purpose. For versions with a lock, the door can be unlocked and opened in an emergency by pressing the lock inwards. After an emergency opening, the lock can simply be reinserted.

### 4.2 Valve inserts

By operating the ball valve, the gas flow can be opened (lever in the direction of flow) or shut off (lever perpendicular to the direction of flow). The pressure or vacuum of the supplied area can be read on the pressure gauge. The NIST connection allows an emergency supply to be established in the event of a repair or a fault in the central gas supply. The emergency inlet point is fitted with a protective cap and a non-return valve, which only opens when a connector is inserted.

#### **Note:** Display of vacuum values

When referring to vacuum values, a distinction is made between absolute pressure and relative pressure (positive or negative pressure relative to ambient pressure). Absolute pressure refers to absolute vacuum (0 bar(a)). Relative pressure, on the other hand, refers to the respective ambient pressure. In the Ventus evo 2, vacuum values are always displayed and monitored as relative pressure. Here, 0 bar corresponds to ambient pressure. A displayed value of -1.0 bar represents, approximately, the maximum achievable vacuum relative to ambient pressure. When setting and assessing limit values for vacuum, it should be noted that all displayed values refer to relative pressure.

### 4.3 Zone Monitor


The range monitor with touchscreen has several menu levels. In the main view, each channel is displayed with its current pressure value. The buttons on the menu bar can be used to access the logbook or enter the settings menu. The menu language is, as per the customer's request, German, English, French, Dutch, Italian, Danish, Latvian or Hungarian.

See also 7. Symbols

#### 4.3.1 Alarm and warning messages

Depending on the configuration, the zone monitor responds to events such as deviations in pressure values or the door being opened with visual and/or audible alerts. An alert remains displayed for as long as the cause of the event persists. Once the cause has been rectified, the alert clears.

'Green' display: Everything is in order. No action is required.


'Yellow' display: An event has occurred that requires action to be taken in the foreseeable future.

'Red' display: An event has occurred that requires immediate action (e.g. line pressure outside the limit values). In the event of an alarm, an acoustic alarm also sounds, which can be silenced by pressing the mute button for up to 15 minutes, depending on the setting.

**Note:** When monitoring vacuum, it may happen that the configured limit (- 1.0 bar) is exceeded and a warning 'Vacuum too high' is displayed. This event is documented in the logbook, but no audible signal is emitted and the warning is not transmitted!

#### 4.3.2 Logbook

The logbook lists all events that have occurred during operation in chronological order. Events include alarms, warnings and changes to device settings. The arrow keys can be used to navigate to older or earlier events. Each event is assigned a sequential number, date and time, and an alarm text. Alarms are displayed in red, warnings in yellow, information in green and system messages in white. For events relating to the analogue or digital inputs, these are assigned to the respective channels via the configuration document, which forms part of the product documentation.


#	Time	Text
558	2021-04-16 13:24	alarms acknowledged
557	2021-04-16 13:24	Added Alarm analog 1 short
556	2021-04-16 13:19	Removed Alarm analog 5 short
555	2021-04-16 13:19	Removed Alarm analog 2 short
554	2021-04-16 13:19	Removed Alarm analog 1 short
553	2021-04-16 13:19	Added Alarm analog 5 short
552	2021-04-16 13:19	Added Alarm analog 2 short

### 4.3.3 Settings

The Settings menu displays the configured parameters for the zone control and shut-off unit, as well as for the individual channels. Changes can only be made using the CAScon USB-RS-485 cable. Use the arrow keys to switch between the individual channel and device settings. Pressing the Home button returns you to the main view.

### 4.3.4 Setting the date and time

The date and time are set in the Settings menu. Use the arrow key to navigate to the last page of the menu, then select the clock icon. Set the date and time using the arrow keys and confirm your entry by pressing the clock icon again. Then return to the main screen.

## 5. Cleaning

The external surfaces of the Ventus evo 2 can be cleaned with a soft, slightly damp cloth. Commercially available cleaning agents and disinfectants may be used for cleaning and disinfection. Only cleaning agents and disinfectants that are suitable for plastic and metal surfaces shall be used.

## 6. Maintenance and repair

Every Ventus evo 2 must undergo a visual inspection at least once a year. During this inspection, the device must be checked for external or obvious damage, the legibility of the display, and the plausibility of the displayed values. In addition, a test run must be carried out by pressing the Test button. The test button is located on the last page under 'Settings'. When the test run is carried out, all the electronic relays are activated, which may trigger alarms in central monitoring systems.

The battery must be replaced every 5 years or in the event of premature discharge. Please follow the instructions in the battery spare parts kit!

Maintenance and repairs must only be carried out by authorised and qualified personnel. Use only genuine Greggersen spare parts for repairs. After repairs, the valve inserts must be checked in accordance with DIN EN ISO 7396-1.

### 6.1 Battery replacement

Loosen the fixing screws on the electronic unit and pull it forwards so that it is still positioned almost vertically. Feel for the battery insertion handle by hand in the battery slot at the bottom right of the battery housing and pull it out.

Feel for the battery slot at the bottom right of the battery housing (the only opening on this side). Insert the replacement battery assembly into the slot with the positive terminal facing the wall until the insertion handle is flush with the electronics housing.

Then swing the electronics unit back into its original position and tighten the fixing screws. Finally, check that the electronics unit and the battery assembly are correctly seated.

### 6.2 Replacing the pressure gauge

Hold the pressure gauge in place, loosen the fixing screw with an Allen key and pull the pressure gauge out. A small amount of gas will now escape through a small hole in the connection. Separate the adapter from the pressure gauge, screw in the new pressure gauge, insert the pressure gauge into the block and retighten the fixing screw. Then carry out a functional and leak test.

### 6.3 Replacing the pressure sensor


Shut off the gas supply by closing the ball valve and relieve the pressure. Disconnect the signal cables from the electronics. Unscrew the sensor from the valve block while counterholding the valve block. Screw the new sensor with a new seal into the valve block while counterholding the valve block. Connect the signal cable to the electronics. Open the gas supply and check the connection for leaks. The Ventus evo 2 can now be put back into operation.

### 6.4 Installation of physical separation

Shut off the gas supply and relieve the pressure. Remove the pressure gauge and loosen the screw connection. Remove the copper gasket, fit the physical isolation device and retighten the screw connection. Use a new copper gasket when reassembling. Once complete, check the connections for leaks. The Ventus evo 2 can now be put back into operation. The physical isolation device is not included with the product but can be purchased separately.


### 6.5 Programming

The Ventus evo 2's electronics can be programmed using the GAScon software and a CAScon USB-RS-485 cable. Please refer to the separate parameterisation instructions.

## 6.6 Aligning the electronics

To adjust the gap between the electronics and the door, first loosen the retaining screws (1). The position of the electronics can now be adjusted using the adjustment screws (2). Then retighten the screws (1).


## 6.7 Disposal

This product contains metal and electronic components. At the end of its service life, it must be disposed of in accordance with the applicable local regulations for waste electrical and electronic equipment. Disposal via household waste is not permitted. Materials should be recycled wherever possible.

## 6.8 Spare parts and accessories

Order No.	Description
902268	Surface-mounted cover
902269	Mounting bracket
905294	CAScon USB-RS-485 cable
905295	CAS CR2032 battery
905292	Basic relay module
905293	Extension relay module
905289	Spare key for Ventus
900497	Ventus repair kit – blind plugs
906595	Ventus 22 physical isolation kit (10 pcs)
905210	Ventus 22 plug-in sensor adapter
905296	Repair set pressure gauge - Ventus 22 O <sub>2</sub>
905297	Repair set pressure gauge - Ventus 22 neutral
905298	Repair set pressure gauge - Ventus 22 Vac
905299	Repair set sensor Ventus 0-16 bar
905300	Repair set sensor Ventus -1 - 6 bar
106042	Schaffner FN2030A-1-06 mains filter
905304	CEGAnet power supply unit 24 V DC/0.5 A
806490	Pipe label gas type, flow direction O <sub>2</sub>
806491	Pipe label gas type, flow direction: AIR
806492	Pipe label gas type, flow direction: VAC
806493	Pipe label gas type, flow direction: N <sub>2</sub> O
806494	Pipe label gas type, flow direction: CO <sub>2</sub>
806495	Pipe label gas type, flow direction: AIR-800
806496	Pipe label gas type, flow direction: AIR <sub>7</sub>
806497	Pipe label gas type, flow direction: N <sub>2</sub>
806498	Pipe label gas type, flow direction: MA
806499	Pipe label gas type, flow direction: SA

## 7. Symbols

	Previous page	Go to the previous page
	Next page	Go to the next page
	Settings	Access system and channel settings
	Test	Automatic system test
	Home	Go to the main view
	Logbook	Access the logbook
	Audio pause	Pause the audible alarm
	Attention	New alarm messages
	Cable break	Cable is broken
	Short circuit	Cable is damaged


## 8. Labelling

	Follow instructions for use		Do not use oil
	Order number		Temperature range
	Medical device		Product contains lead
	Manufacturer		

### Labels for Ventus evo 2

<b>REF 902252</b> Ventus evo 2, 4-6fach mit Fenster und Schloss Ventus evo 2, 4-6fold with window and lock  (01) 04031196020493 (11) 260423 (21) 902252/00001 Greggersen Gasetechnik GmbH Bodestr. 27-29 21031 Hamburg Germany	 <small>www.greggersen.de/identifizierte</small> 	Teil der Konfigurationsbaugruppe Part of the configuration assembly <b>Ventus evo 2, 4-6fach FS</b> SN 902252/00001 Greggersen Gasetechnik GmbH
<b>REF 902252</b> Ventus evo 2, 6-fach Input: 230V/50Hz, <24W  (01) 04031196020493 (11) 260423 (21) 902252/00001 Greggersen Gasetechnik GmbH Bodestr. 27-29 21031 Hamburg Germany	  CE 0123	

### Labels for Ventus valve inserts

<b>REF 900780</b> Ventus 22 Ventileinsatz valve insert for Ventus 22 medium: O2  (01) 040311960012269 (11) 260129 (10) 12345 Greggersen Gasetechnik GmbH Bodestr. 27-29 21031 Hamburg Germany	 <small>www.greggersen.de/identifizierte</small> 	<b>REF 900780</b> Ventus 22 Ventileinsatz medium: O2 (10) 12345 Greggersen Gasetechnik GmbH
--	---	---

## 9. Technical data

Classification according to (EU) 2017/745:	active, non-invasive, Class IIa
Equipment protection class:	1
Protection rating:	IP 40
Configuration:	up to six valve inserts comprising ball valve, pressure gauge, pressure sensor, NIST connection
Monitoring:	Range monitor in accordance with DIN EN ISO 7396-1
Signal outputs:	Optional relay modules with potential-free contacts CEGANET ModBus module Ethernet interface
Supply voltage:	230 V AC / 50 Hz
Power supply:	24 V DC, 0.625 A
Switching capacity of relay module:	30 V DC, 1 A or 125 V AC, 0.6 A
Piping:	Cu pipe 22 x 1 mm (top feed)
Nominal distribution network pressure:	400...500 kPa (compressed gases) 700...1,000 kPa (pneumatic tools) < - 40 kPa (vacuum)
Test pressure: max.:	1,600 kPa
Weight:	3-port: 7.5 kg, 6-port: 11 kg
Dimensions:	see page 182
Operating temperature:	+10 °C to +40 °C
Storage temperature:	-10 °C to +70 °C
Relative humidity:	15 % – 95 %
Atmospheric pressure:	700 hPa – 1,060 hPa
Expected service life:	10 years

## 10. Error messages

<b>Error message</b>	<b>Possible cause</b>	<b>Action</b>
Battery flat	Battery has run flat or it is past its best	Replace the battery
No network connection	Ethernet cable has no	Check connection
Operating voltage too low	The power source has less than 22 V (e.g.: power cut, power supply unit faulty)	Check power supply
Internal error	Incorrect value set in the configuration Error in the programme code Base module	Notify service
RM bus error	Error in the programme code Relay module	Notify service
Cable break	Connection cable disconnected Sensor faulty Incorrect value set	Check connection Replace sensor Notify service
Short circuit	Short circuit Connection cable Sensor faulty Incorrect value set	Check connection Replace sensor Notify service
RM communication error	Connection cable interrupted	Check connection

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