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

**Dear customers,**

we would like to thank you for purchasing this Greggersen Gasetechnik GmbH product. For questions and information, our sales and support team will be happy to assist you.

+49-(0)40 739 357-0; [sales@greggersen.de](mailto:sales@greggersen.de)

Please read the instructions for use carefully before using the medical device. Use it only if you have fully understood the application and function. All serious incidents must be reported to Greggersen Gasetechnik GmbH as well as to the health authority responsible for the place of residence.

**Medical Device in accordance with:**

DIN EN ISO 7396-2

DIN EN ISO 9170-2

MDR 2017/745

**Foreword**

The AGSS Forano terminal unit of Greggersen Gasetechnik GmbH complies with the relevant provisions of Regulation (EU) 2017/745 of the European Parliament and of the Council (MDR 2017/745) and is classified as an active, non-invasive Class I medical device.

**Intended Purpose**

For the removal and transfer of excess residual gases during inhalation anesthesia to a suitable disposal system.

**Principle of Operation**

Terminal units in Anaesthetic Gas Scavenging Systems (AGSS) are connected directly to a central compressed gas supply with a supply pressure of approximately 500 kPa compressed air or nitrogen.

The ejector is connected to the gas supply and generates negative pressure during operation. This negative pressure extracts excess anaesthetic gases and conveys them via the suction line together with the compressed gas to the outside atmosphere or into the exhaust duct of an air-conditioning system.

The terminal unit is activated by connecting (coupling) an AGSS plug or a hose with a corresponding connector. When the connector is disconnected, the suction process stops.

A visual indicator is integrated for operational monitoring.


## Warnings

- Prior to using the device, the user must be familiar with the product and the instructions for use and must have understood them.
- Before each use, the terminal unit must be inspected for external damage or contamination. It must be ensured that the terminal unit is correctly aligned and that the labeling is positioned horizontally.
- During cleaning, ensure that no liquids enter the terminal unit. After cleaning, the terminal unit must be completely dried.
- After installation of the terminal unit, an acceptance test in accordance with DIN EN ISO 7396-2 must be carried out.
- Only use the terminal unit with AGSS connectors of type 1L.
- All serious incidents must be reported to Greggersen and to the competent authority.

## Drawing and Parts List

1. O-ring 2.5 x 1.3
2. AGSS closing cone
3. O-ring 9.25 x 1.78
4. O-ring 14 x 2
5. Sealing washer 25 x 18 x 1
6. O-ring 16.5 x 1
7. Forano pressure ring AGSS
8. Driving gas connection
9. Visual indicator connection
10. Ejector driving gas connection

## Accessories

- 902071 AGSS plug
- 902095 AGSS plug, 45° angled

## Patient Target Group

The patient has no direct contact with the product. Therefore, there are no restrictions regarding the patient target group.

When the AGSS Forano terminal unit is used as intended, the focus is not on patient protection but on protecting persons in the immediate surroundings, in particular medical staff.

## Intended Users

The medical device is intended for clinical use by medically trained healthcare professionals (therapists, nursing staff, physicians).

## Operation

Before each commissioning, a visual inspection and functional test must be carried out. Following visual inspection for obvious damage or contamination, the pressure ring must be actuated to check smooth operation.

To establish a connection, an AGSS plug or a hose with a type 1L AGSS connector is inserted into the terminal unit.

By pressing the pressure ring of the terminal unit, the plug is released from the terminal unit. The plug should be firmly held by hand during this process.

Suction is interrupted.

### **Cleaning**

The terminal unit may be cleaned with commercially available disinfectants and cleaning agents using a damp method. Ensure that no liquids enter the terminal unit during cleaning. Before reuse, the terminal unit must be completely dry. The hygiene plan of the healthcare facility must be observed.

### **Installation**

- The installation housing (cavity wall / flush-mounted version) must first be secured in the wall.
- Cavity wall mounting: Secure the housing with 6 screws on the left and right sides; solder pipes to the pipeline system.
- Flush-mounted installation: Plaster the housing into the wall; solder pipes to the pipeline system.
- Connect the terminal unit to the pipe or ejector connections using the connecting lines (the ejector is plug-in type and secured with a nut).
- Insert the terminal unit with the front plate into the housing and secure with 2 screws.
- Supply systems: Secure the terminal unit with 2 screws, insert the ejector into the exhaust line and secure with a nut.
- Connect the terminal unit, ejector, and pipe connections using the connecting lines.
- Adjust suction performance using the valve spindle (11 o'clock position) with a 2.5 mm hex key.
- Press the pressure ring onto the terminal unit.
- Test the Anaesthetic Gas Scavenging System in accordance with DIN EN ISO 7396-2.

### **Maintenance and Inspection**

Terminal units must undergo visual inspection (external damage) before each use. At least once per year, a functional and leak test must be carried out (replace wear parts if necessary). Every six years, maintenance must be performed (metrological and safety inspection). In areas with particularly frequent use, maintenance intervals must be shortened accordingly.

Malfunctions and damage may only be rectified by authorized qualified personnel.

Only original GREGGERSEN spare parts may be used for repairs.

### **Disposal**

The terminal unit must be disposed of in accordance with the applicable national regulations.

## Technical Data


Design:	All-metal construction, plug-in ejector for vacuum generation, one-handed operation for coupling and uncoupling, visual indicator for operation control, integrated valve for performance adjustment.
Connections:	<u>Supply lines for flush-mounted / cavity wall / surface-mounted systems</u> Compressed gas line: Cu tube 8 x 1 Exhaust line: Cu tube 15 x 1 <u>Terminal unit / Ejector</u> Compressed gas supply: Coupling for 4 mm PA hose or 7 mm hose nozzle Ejector: Coupling for 4 mm PA hose or 7 mm hose nozzle Anaesthetic gas outlet: 15 mm spiral hose Visual indicator: Coupling for 4 mm PA hose
Performance:	Minimum 50 l/min at 500 kPa compressed gas
Gas consumption:	Approx. 12 l/min at a suction performance of 50 l/min
Installation:	Flush-mounted or cavity wall installation, Ceiling-mounted supply systems, Wall-mounted supply systems


Our products are partly manufactured from brass containing more than 0.1% lead. Lead is metallurgically bound within the alloy, so no exposure is expected.

## Installation Dimensions


### Flush-mounted / cavity wall


### Surface-mounted


## Rail system

## Ceiling-mounted supply systems (DVE)

## Visual Indicator

## Maintenance Parts

Article Number	Description
900442	Repair kit AGSS Forano
900058	Ejector
100058	Visual indicator
805858	Pressure ring Forano AGSS
100055	Hose connector 4 mm 90° or hose nozzle 6 mm for driving gas connection, visual indicator connection and ejector driving gas connection

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