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
# Ventus evo

## Dear customers,

Thank you for purchasing this Greggersen product.  
If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739357-0, [sales@greggersen.de](mailto:sales@greggersen.de)

Please read the instructions for use thoroughly before using the medical device. Only use it once you have fully understood the application and function!

## 1. General information

The Ventus evo product from Greggersen Gasetechnik GmbH complies with the requirements of Regulation 93/42/EWG and is classified as an active and non-invasive Class IIb device.

### 1.1 Medical device in accordance with

DIN EN ISO 7396-1

DIN EN 60601-1

93/42/EWG

### 1.2 Purpose

For monitoring and displaying the pressure data of various medical gases or vacuum of a ward or section, for logging warnings and alarms, for disconnecting the gas or vacuum network of a zone in the event of emergency or maintenance, and for emergency supply of a zone.

### 1.3 Intended users

The intended users are medically or technically trained staff.

### 1.4 Patient target group

All patient groups that need to be supplied with medical gases or vacuum or need to be treated with gas-operated treatment devices. The product is used to provide gases or vacuum. Its use is decided by the medical staff.

### 1.5 Warnings

- Before using the device, the user must have understood the instructions for use and needs to be familiar with the product.
- The device may only be used by persons who have been briefed and are familiar with the operating procedures of the respective healthcare facility.
- All serious incidents must be reported to Greggersen Gasetechnik GmbH and to the health authority responsible for the place of residence.

## 2. Device description


The respective zone can be separated from the medical gas supply network using a valve insert.

The Ventus evo zone shut-off units are configured according to the customers requirements for the respective zone.

All zone shut-off units consist of a housing made of powder-coated metal, which can be opened via a door. The door is available with and without a viewing window, as well as with and without a lock. By pushing in the lock or the disc, it is possible to open the door, in case of an emergency. Zone cut-off units without a lock can always be opened by pulling on the right side of the door.

The system is available in two sizes for 1-3 or 4-6 media and alternatively with a surface-mount frame.

According to the requirements of the ward, the zone shut-off units are equipped with valve inserts for the different types of gas. The valve inserts allow the gases to be passed one way from the source to the zone. A ball valve can be used to separate the zone from the source at this point. The valve inserts are equipped with a NIST connection behind the ball valve and optionally in front of the ball valve. Emergency care can be provided from this emergency intake point. A pressure sensor is mounted in the valve insert behind the ball valve, which continuously measures the gas pressure of the respective zone and transmits it to the zone monitor.

The measurement data are processed in the zone monitor and compared with the individually specified limit values of the health facility. In the event of a fault, the staff are alerted visually and acoustically. Optionally, the door can be monitored via the zone monitor and the signal can be made available via a potential-free contact. To forward the various signals and warnings to a building control system, the main board can be expanded with up to five relay modules.

In the event of a power failure, this is indicated by a signal lamp on the front panel.

### 3. Installation and starting operation

Installation and starting operation must be carried out by qualified personnel.

Before starting operation for the first time or resuming operation, the product must be tested in accordance with the respective national regulations for electrical equipment.

#### 3.1 Checking delivery

Prior to installation, the product must be checked for damage to the transport packaging, for obvious damage to the product, and for the completeness and correctness of the configuration according to the delivery note.

#### 3.2 Required tools and aids

To install the Ventus evo zone shut-off unit, only ordinary tools, such as drills, cordless screwdrivers, Allen keys, spirit levels, wire strippers and multimeters, are required. A spatula and filler compound may be required to close the gaps.

#### 3.3 Installation

Remove the housing from the packaging and disassemble the electronics. To do this, remove the monitor cover by removing the fastening screws. Then loosen the screws above the display, tilt the electronics forward and remove by lifting slightly. Cut free the bushings for electrical cables and for the copper pipes in the rear section.

##### 3.3.1 Concealed installation in concrete or masonry walls

Chisel out the installation opening for the housing and install the housing horizontally in the wall through the rear or side holes. The front lower edge of the housing must be flush with the modular wall. Plaster layer thickness, tile level, etc. must be taken into account here.

Install the valve assemblies in rear section according to the configuration.

Insert the electrical cables into the housing.

##### 3.3.2 Installation in plasterboard walls

Establish the installation opening in the plasterboard wall. The opening should be no more than 5 mm larger than the housing on all sides. Where appropriate, take fire protection enclosures into consideration.


Fix the housing horizontally to the stud frame in the wall using the holes in the housing.

Insert the electrical cables into the housing.

##### 3.3.3 Surface-mount installation

Install the surface-mount frame horizontally on the wall, push in the rear section and screw tight. Install the valve assemblies in rear section according to the configuration.

Insert the electrical cables into the housing.

### 3.4 Gas connection

During the soldering process, open the ball valve and purge the pipes with inert gas. The effect of heat must be taken into account and, if necessary, the housing and pipes must be protected. The rear section is then covered with the plaster protection. Before closing the walls, a pressure test of the system in accordance with DIN EN 7396-1 must be carried out. For vacuum valve inserts, it should be noted that the pressure test must be carried out without a manometer, otherwise it would be damaged. The valve block is delivered fitted with a dummy plug for this purpose. Following the pressure test, remove the dummy plug and insert the manometer. See 6.1 Replacing the manometer.


#### Apply labels to indicate flow direction

After soldering the pipes, the labels to indicate the direction of flow and the type of gas must be glued to the left-hand leg of the valve blocks. For compressed gases, the label must be attached so that the arrow points upwards; for vacuum, it must point downwards.

### 3.5 Electrical connection

The electrical connection may only be carried out in a de-energised state! Strip the cable ends 8-9 mm. Use wire-end ferrules for flexible cables.


#### Suggestion for cable selection


Mains voltage terminal blocks:	NYM- 3x1.5 mm <sup>2</sup> - bis 3x2.5 mm <sup>2</sup>
Analog/Digital inputs:	LiCY- 2x2x 0.5 mm <sup>2</sup> (flexible) Do not connect shield IStY- 2x2x0.6 mm <sup>2</sup> (rigid)
Relay modules:	IStY- 2x2x0.6 mm <sup>2</sup> IStY- 6x2x0.6 mm <sup>2</sup>
Network cable:	at least CAT 6

#### 3.5.1 Installation of the electrical components

1. Install the terminal blocks (1) on the DIN rail (2)
2. Fasten mains filter (3) with screws
3. Clamp power supply unit (4) to the DIN rail
4. Install connecting cables (5)

**Note:** By plugging in the mains plug and closing the jumper, the optical alarm is enabled with an LED in the event of mains voltage failure. To prevent the internal battery from discharging too quickly, disconnect the mains plug or remove the jumper during transport, while starting operation or other prolonged periods without power supply.

### 3.5.2 Connection of sensors and relay module

1. Connect the sensor cables to the electronics. The terminal assignment is shown on the back of the electronics.
2. Fasten the mains filter with screws.
3. Clamp the power supply unit to the DIN rail.
4. Install the connecting cables.


### 3.5.3 Terminal assignment of the main board

Pin	Function
1	Channel 4
2	Channel 4
3	Channel 3
4	Channel 3
5	Channel 2
6	Channel 2
7	Channel 1
8	Channel 1

Pin	Function
1	Channel 8
2	Channel 8
3	Channel 7
4	Channel 7
5	Channel 6
6	Channel 6
7	Channel 5
8	Channel 5

Pin	Function
1	Channel 3 -
2	Channel 3 +
3	Channel 2 -
4	Channel 2 +
5	Channel 1 -
6	Channel 1 +

Pin	Function
1	Channel 6 -
2	Channel 6 +
3	Channel 5 -
4	Channel 5 +
5	Channel 4 -
6	Channel 4 +

Pin	Function	Remark
1	+ 24V	
2	GND	
3	Jumper battery	Please refer 3-5.1
4		

Pin	Function	Color cable
1	24 VDC	red
2	D-	white
3	D+	green
4	GND	black


Pin	Function
1	NC
2	C
3	NO

### 3.5.4 Terminal assignment of the relay module


Relay 1-6

Pin	Function
1	Relay 3
2	Relay 3
3	Relay 2
4	Relay 2
5	Relay 1
6	Relay 1

Relay 1-3


Pin	Function
6	Relay 4
5	Relay 4
4	Relay 5
3	Relay 5
2	Relay 6
1	Relay 6

Connection relay module main board

#### Setting the bus addresses

RM	S <sub>1</sub>	S <sub>2</sub>	S <sub>3</sub>	S <sub>4</sub>
1	OFF	OFF	OFF	OFF
2	ON	OFF	OFF	OFF
3	OFF	ON	OFF	OFF
4	ON	ON	OFF	OFF
5	OFF	OFF	ON	OFF

DIP switch

## 4. Operation

The Ventus evo zone shut-off unit is designed to be operated intuitively and by untrained personnel in an emergency.

#### 4.1 Opening the door

In an emergency, such as a gas leak or fire, or for maintenance, it may be necessary to open the door to access the inside of the shut-off unit. The version without lock can be opened without further ado, while the version with lock can be opened with the appropriate key. To open the door, pull it open on the right and open it fully. The door can also be opened in an emergency by pushing in the lock or the disc. The disc and lock are easy to re-insert after emergency opening.


#### 4.2 Valve insert

By actuating the ball valve, the gas flow can be opened (lever in the flow direction) or interrupted (lever against the flow direction). The pressure or vacuum of the zone supplied can be read off the manometer. The NIST connection allows an emergency supply in case of repair or failure of the central gas supply. The emergency inlet point is fitted with a protective cap and a non-return valve which only opens when a plug is inserted.

#### 4.3 Zone monitor


The zone monitor with touchscreen has several menu levels. In the main view, each channel is displayed with the current pressure value. The keys in the menu line can be used to access the **Logbook** or to enter the **Settings range**. According to the customer's requirements, the menu language can be German or English.

Also see 7. **Symbols**

##### 4.3.1 Alarm and warning signals

Depending on the configuration, the monitor responds to events such as deviating pressure values or opening of the door, with visual and/or acoustic signals. A signal is displayed as long as the cause of the event continues. If the cause is eliminated, this signal also disappears.

Display „green“: Everything is OK. No actions are necessary.

Display „yellow“: An event has occurred that requires necessary action in the foreseeable future.

Display „red“: An event has occurred that requires immediate action (e.g. line pressure outside the limit values). In the event of alarms, an acoustic alarm also sounds, which, depending on the setting, can be suppressed for up to 15 minutes by pressing the mute button.


**Note:** When monitoring vacuum, it may occur that the lower limit value (- 1.0 bar) is exceeded and a warning „Vacuum too high“ is displayed. This event is documented in the logbook, but there is no acoustic signal tone and this warning is not forwarded!

### 4.3.2 Logbook

The logbook chronologically lists all events that occurred during operation. Events are alarms, warnings or also changes of device settings. Use the arrow keys to navigate to older or earlier events.

Each event is assigned a consecutive number, date and time, as well as an alarm text.


Alarms are displayed in red, warnings in yellow, informations in green, and system signals in white. For analog or digital input events, an assignment to the respective channels is made with the configuration document, which is part of the product documentation.


#	Time	Text
558	2021-04-16 13:24	alarms acknowledged
557	2021-04-16 13:24	Added Alarm analog 1 short
556	2021-04-16 13:19	Removed Alarm analog 5 short
555	2021-04-16 13:19	Removed Alarm analog 2 short
554	2021-04-16 13:19	Removed Alarm analog 1 short
553	2021-04-16 13:19	Added Alarm analog 5 short
552	2021-04-16 13:19	Added Alarm analog 2 short

**4.3.3 Settings** The set parameters of the monitor and the individual channels are displayed in the Settings menu. Changes are only possible with the CAScon cable USB-RS-485.

Use the arrow keys to switch between the individual channel or device settings. Press the Home button to switch back to the main view.

SETTINGS		8 of 8 Ventus
DEVICE		
Next service	2026-03-27	
MAC adress	04:91:62:B6:D1:1C	

SETTINGS		6 of 8 Ventus
CHANNEL 6		
Name	Door	Type
Low value		Normal state
High value		Alarm state

### 4.3.4 Set date and clock

The clock-setting is located under SETTINGS. With arrow key (right arrow) to last page. Click on the clock here. Press arrow up or down to set the clock. When the setting is complete press against clock symbol as confirmation. Back to main screen.

## 5. Cleaning

The outside of the Ventus evo can be wiped with a slightly moistened soft cloth. Commercially available cleaning agents and disinfectants can be used.

Do not use cleaning agents or disinfectants that are not suitable for common plastic or metal surfaces.


## 6. Maintenance and repair

Every Ventus evo must be visually inspected at least once a year. The device must be checked for external or obvious damage, legibility of the displays and plausibility of the displayed values. In addition, a test run must be carried out by pressing the Test button. The Test button is located on the last page under Settings. During the test run, all relays of the electronics are enabled, which can trigger alarms for central monitoring.

The battery must be replaced every 5 years or in the event of premature discharge. The information in the battery spare parts kit must be observed for this purpose!

Maintenance and repairs must only be carried out by authorised expert personnel. Only use genuine Greggersen spare parts when performing repairs. After repairs, the valve inserts must be checked in accordance with DIN EN ISO 7396-1.

Pos.	Description
1	Inlet pipe
2	Copper gasket
3	Ball valve G3/4"
4	Physical separation
5	Fastening hole for installation
6	Pressure sensor
7	NIST emergency intake point
8	Pressure gauge fastening screw
9	Gauge
10	Branch pipe

### 6.1 Replacing the manometer

Hold the manometer firmly, loosen the fastening screw with Allen key and pull out manometer.

Some gas now escapes through a small hole in the connection.

Separate the adapter and manometer, screw in a new manometer

Insert the manometer into the block and retighten the fastening screw. Then carry out a function and sealing test.

### 6.2 Replacing the pressure sensor

Shut off the gas supply by closing the ball valve and release the pressure. Disconnect the signal lines from the electronics. Screw the sensor out of the valve block and lock the valve block. Screw the new sensor with new gasket into the valve block while countering the valve block. Connect signal line to the electronics. Open the gas supply and check the connection for leaks. The Ventus evo can be put back into operation.

### 6.3 Installation of physical separation


Interrupt the gas supply and relieve the pressure, remove the manometer and loosen the screw connection. Remove the copper gasket and install physical separation then retighten the screw connection. Use a new copper gasket after dismantling. After completion, check connections for leaks. The Ventus evo can be put back into operation. The physical separation is not part of the product, but can be purchased separately.

### 6.4 Programming

The electronics of the Ventus evo can be programmed using of GAScon software and the CAScon cable USB-RS-485. The separate parameterization instructions must be observed here.

### 6.5 Aligning the electronics

To adjust the gap between the electronics and the door, first loosen the retaining screws (1). The position of the electronics can now be changed with the adjusting screws (2). Then the screws (1) must be tightened again.

### 6.5 Spare parts and accessories

Item No.	Description
905294	CAScon USB-RS-485 cable
905295	CAS CR2032 battery
905292	Relay module, basic
905293	Relay module, expansion
905289	Replacement key for Ventus
900497	Ventus repair kit, blind plug
906595	Physical separation Ventus 22, 10 pcs.
905296	Repair kit Mano Ventus 22 O <sub>2</sub>
905297	Repair kit Mano Ventus 22, neutral
905298	Repair kit Mano Ventus 22 VAC
905299	Repair kit Sensor Ventus 0–16 bar
905300	Repair kit Sensor Ventus -1 to 6 bar
106042	Mains filter Schaffner FN2030A-1-06
905304	CEGANet power supply 24V DC / 0.5A

## 7 Sybols

	Previous page	Go to previous page
	Next page	Go to next page Ventus
	Settings	Calls up the system and channel settings
	Test	Automatic system test
	Home	Go to main view
	Logbook	Call logbook
	Audio pause	Interruption of the acoustic alarm
	Caution	new alarm signals
	Cable break	Cable is broken
	Short circuit	Cable is damaged

## 8. Technical data

Device protection class:	1
Basis	UDI:4031196VENTUSUH
Protection type:	IP 50
Configuration:	up to six valve inserts consisting of ball valve, manometer, pressure sensor, physical separation, NIST connection
Monitoring:	Zone monitor in accordance with DIN EN ISO 7396-1
Outputs:	potential-free relay (collective signal), Ethernet interface
Optional:	relay module (up to 30 potential-free contacts)
External supply voltage:	110 VAC - 230 VAC (50 - 60 Hz)
Internal supply voltage:	24 V DC
Power consumption electronics:	max. 0.5 A
Switching capacity relay module:	30 V DC / 1 A or 125 V AC / 0.6 A
Pipe:	Cu pipe 22x1 (infeed from above)
Nominal operating pressure:	400...500 kPa (compressed gas) 700...1,000 kPa (pneumatic tools) ◁ - 40 kPa (vacuum)
Test pressure:	max.: 1600 kPa

### Concealed / cavity wall housing

Installation housing:	1-3x (HxWxD) approx: 470x328x69 mm
Ventus evo cover:	1-3x (HxWxD) approx: 517x390x17 mm
Installation housing:	4-6x (HxWxD) approx: 470x580x69 mm
Ventus evo cover:	4-6x (HxWxD) approx: 517x630x17 mm

### Surface-mounted housing

1-3 housings (HxWxD) approx:	517x390x102 mm
4-6 housings (HxWxD) approx:	517x630x102 mm
Operating temperature:	+10°C...+40°C
Storage temperature:	-10°C...+70°C
Humidity:	15%...95%
Air pressure:	700 hPa...1,060 hPa

## 9. Error messages

Error messages	Possible causes	Action
Low battery	Battery is empty or too old	Change battery
No network	Ethernet line has no connection	Check connection
Low supply	The voltage source is less than 22 V (e.g.: power failure, power supply defective)	Check power supply
	Incorrect value set in the configuration	Notify Service
Internal error	Error in program code basic module	Notify Service
Bus error	Error in program code relay module	Notify Service
Cable cut	Connection line interrupted	Check connection
	Sensor defective	Change sensor
	Incorrect value set	Notify Service
Short circuit	Short circuit connecting cable	Check connection
	Sensor defective	Change sensor
	Incorrect value set	Notify Service
RM comm. error	Connecting cable interrupted	Check connection

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