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

# Flowmeter Fink

**Dear customers,**

we would like to thank you for purchasing this Greggersen Gasetechnik GmbH product. For questions and information, our sales and support team will be happy to assist you.

+49-(0)40 739 357-0; [sales@greggersen.de](mailto:sales@greggersen.de)

Please read the instructions for use carefully before using the medical device. Use it only if you have fully understood the application and function. All serious incidents must be reported to Greggersen Gasetechnik GmbH as well as to the health authority responsible for the place of residence.

**Medical device according to:**

DIN EN ISO 15002  
DIN EN ISO 19054  
(EU) 2017/745

**Connections according to:**

DIN 13260-2	SS 875 24 30 (AGA)	Medap
BS 5682	SN ENV 737-6:2005	NF S90-116
DIN EN ISO 18082		

**Symbols**

Observe instructions for use

Do not use oil

Date of manufacture, calendar week – year

Product contains lead


Manufacturer

Order number

**Labels**

<b>REF 900753</b>	eIFU 
Flowmeter, DIN, standard	<a href="http://www.greggersen.de/downloads/">www.greggersen.de/downloads/</a>
Flowmeter, DIN, standard	
O2, 0-15 l/min, Stecker / PLUG-UNIT	
 (01) 04031196001508	<b>Lot</b> 826
(11) 260217	<b>CE 0482</b>
(10) LOT826	
 Greggersen Gasetechnik GmbH Bodestr. 27-29 21031 Hamburg Germany	

<b>Medical Flowmeter</b>	
input: 400 ... 500 kPa	
output: 0 ... 15 l/min	<b>O<sub>2</sub></b>
SN900753/20104	<b>CE 0482</b>
<b>Greggersen Gasetechnik GmbH, 21031 Hamburg</b>	

## Intended purpose

The Fink device is intended for the supply of medical gases with a selected flow rate to a downstream medical device. Through the function of the device, for example, oxygen can be administered to a patient via a face mask or nasal cannula. Plug-in devices are directly connected to a terminal unit; the version with NIST hose connection is mounted on an equipment rail and supplied with pressurized gas via a low-pressure hose. Using an adjustment valve, the flow rate is set to a desired value, which can be read on the ball of the variable-area flowmeter. The medium is provided at an outlet connection. The use in the professional healthcare sector in medical facilities is decided by medical professionals.

## Warnings

- The device may only be used according to its intended purpose.
- Before each use, it must be ensured that the product is correctly aligned.
- The products partly consist of brass with a lead content of more than 0,1%. The lead is firmly bound in the material, therefore no exposure is expected.
- The product is not suitable for sterilization.
- No structural modifications may be made to the product.
- No items that do not serve immediate use may be attached to the product.
- Connect or disconnect NIST connections and accessories only when the system is depressurized.
- Oxygen is a strongly oxidizing gas. Always keep all connections free of oil and grease! The use of hand cream should be avoided during handling. Do not smoke and keep away from open flames near oxygen equipment! Fire hazard!
- All serious incidents must be reported to Greggersen Gasetechnik GmbH as well as to the competent authorities responsible for the place of residence.

## Variants

The Fink flowmeters are available in the following configurations:

- Devices for the gases O<sub>2</sub> (0–15 l/min), AIR (0–15 l/min) and CO<sub>2</sub> (0–12 l/min)
- Plug connections according to standards: DIN, NF, BS, SS, Medap, SN
- Rail-mounted devices with NIST connection; rail clamp DIN and CH
- Single flowmeter and double flowmeter
- Various lengths (distance between device and terminal unit)
- Other versions on request

## Operation


Before each use, the product must be checked for proper function, smooth operation of the valves, as well as for external damage, obvious leaks, and contamination. Non-functional, damaged, leaking, or contaminated devices must not be put into operation. Plug-in devices must be connected directly to a terminal unit. Rail-mounted devices must be fixed to an equipment rail according to DIN EN ISO 19054 before operation, and the low-pressure hose must be connected via the gas-type-coded NIST fitting. Set the device to “0” before connection.

Set the device to the desired value using the handwheel. Check whether the device produces a gas flow that can be heard, felt, or, when using a humidifier, seen. This ensures that an actual gas flow is generated and that the device is not, for example, in the parking position.

Afterwards, a connection to the patient can be made via a hose nipple at the 9/16-18 UNF connection or alternatively via a humidifier. The flow can be set to different values using the handwheel, which can be read on the ball in the measuring tube.

After use, set the flowmeter to “0”. The device may remain in the terminal unit or be moved to the parking position. If necessary, disconnect the humidifier or hose connection from the flowmeter. The device can now be disconnected from the terminal unit and, if applicable, removed from the equipment rail. When using humidifiers, ensure that no liquid from the humidifier enters the device after it has been removed from the terminal unit.

### Cleaning

It is recommended to clean the device after each use by wipe disinfection according to the manufacturer’s instructions, ensuring compatibility with polycarbonate. Commercially available disinfectants should be used. During cleaning, ensure that no liquids enter the product. Before reuse, the device must be completely dry. Do not soak the device with disinfectants! The device may be damaged. The hygiene plan of the facility must be observed. The device is intended for multiple use with different patients.

### Disposal

There are no special disposal requirements for this device. The device may be disposed of with household waste. If the device is contaminated, appropriate precautions must be taken.

### Maintenance and inspection

The device must be serviced at least every 5 years. An inspection including functional and leak testing must be carried out at least every 2 years. Shorter maintenance and inspection intervals may be necessary in case of unusually heavy use. Maintenance and repairs may only be carried out by authorized qualified personnel. Use only original Greggensen spare parts.

### Technical data


Operating pressure:	400–500 kPa
Adjustment accuracy:	±10% of the set value
Outlet:	9/16-18 UNF male thread
Weight:	Single device plug 395 g Double device plug 965 g Single device rail 760 g Double device rail 1690 g
Storage conditions:	–20 ... +60 °C 15 ... 93% relative humidity 620 ... 93% hPa air pressure
Operating conditions:	10 ... +40 °C
Expected service life:	10 years

### Spare parts and accessories

- 900430 repair set flowmeter tube
- 900619 hose connection 9/16-18 UNF
- 904836 humidifier unit 0.25 l
- Low-pressure hoses

### Components

Pos	Description
1	Handwheel
2	Flow tube with ball
3	Gas-specific connector
4	Rail clamp
5	Pressure inlet (NIST)
6	Flow outlet 9/16-18 UNF
7	Hose connection 9/16-18 UNF

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